



MEDICAL

---

## Guide Book

Effective February 1, 2019

## General Information

Ordering Notes.....	4
Freight and Shipping.....	4
Midmark and Ritter Medical Equipment Warranty.....	5

## Midmark Chairs and Tables

Midmark 641 Procedure Chair.....	6-7
Midmark 630 HUMANFORM® Procedure Chair.....	8-10
Midmark 647 Podiatry Chair.....	11-12
Midmark 646 Podiatry Chair.....	13
Midmark 626 Barrier-Free® Examination Chair.....	14-17
Midmark 640 Pediatric Examination Table.....	18
Upholstery Color Palette.....	19

## Ritter Chairs and Tables

Ritter 230 Procedure Chair.....	20-21
Ritter 244 Barrier-Free® Examination Table.....	22
Ritter 225 Barrier-Free® Examination Chair.....	23
Ritter 224 Barrier-Free® Examination Chair.....	24-25
Ritter 204 Manual Examination Table.....	26-27
Ritter 203 Manual Examination Table.....	28
Ritter 95 Manual Examination Table.....	28

## Seating

Midmark 425 Physician Seating (Hand Operated).....	29
Midmark 427 Physician Seating (Foot Operated).....	29
Ritter 279 Physician Seating (Hand Operated).....	30
Ritter 278 Physician Seating (Foot Operated).....	30
Ritter 277 Physician Seating with Backrest.....	30
Ritter 276 Physician Seating.....	30
Ritter 275 Physician Seating with Backrest.....	31
Ritter 274 Physician Seating.....	31
Ritter 273 Physician Seating with Backrest.....	31
Ritter 272 Physician Seating.....	31
Ritter 271 Physician Seating with Backrest.....	32
Ritter 270 Physician Seating.....	32
Ritter 280 Patient Side Chair.....	33
Ritter 282 Patient Bariatric Side Chair.....	33
Midmark 680 Patient Side Chair.....	34
Ritter 281 Blood Draw Chair.....	34
Seating Upholstery Color Selections.....	35

## Lighting

Ritter 255 LED Procedure Light.....	36
Ritter 253 LED Exam Light.....	37
Ritter 250 LED Exam Light.....	38

## Diagnostic Products

Midmark Diagnostic Medical Products Warranty.....	39
Midmark Diagnostic Medical Products Ordering Information.....	40
Midmark HEINE Diagnostic Instruments Limited Warranty.....	41
Midmark Diagnostic Medical Products Support Information.....	42
Midmark Digital Vital Signs and Digital Vital Signs Monitor and Supplies.....	43-47
Midmark Digital Spirometer, Software and Supplies.....	48
Midmark Digital ECG, Software and Supplies.....	49
Midmark Digital Holter, Software and Supplies.....	50-51
Midmark HEINE Diagnostic Instrument Wall Kits and Accessories.....	52
Midmark HEINE Diagnostic Instrument Desktop Kits and Accessories.....	53
Midmark HEINE Diagnostic Instrument Replacement Parts and Supplies.....	54
Midmark HEINE Gamma® XXL LF Sphygmomanometer Kit and Accessories.....	55
Midmark HEINE Gamma® G5 Sphygmomanometer Kit and Accessories.....	56
Midmark Diagnostic Cart Accessories and Supplies.....	57

**Midmark Workstations**

Midmark Workstation Warranty..... 58  
 Midmark 6211 Workstation ..... 59  
 Midmark 6212 Workstation ..... 60  
 Midmark 6213 Workstation ..... 61  
 Midmark 6214 Workstation ..... 62  
 Midmark 6215 Workstation ..... 63  
 Midmark 6217 Workstation ..... 64  
 Midmark 6218 Workstation ..... 65  
 Midmark 6219 Workstation ..... 66  
 Midmark 6281 Workstation ..... 67  
 Midmark 6282 Workstation ..... 68  
 Midmark 6283 Workstation ..... 69

**Midmark Synthesis® Cabinetry Collection**

Midmark Synthesis Cabinetry Warranty..... 70  
 Midmark Synthesis Cabinetry Color Selections ..... 71  
 Midmark Synthesis Exam Room Packages E1-E7 ..... 72-78  
 Midmark Synthesis Mobile Treatment Cabinets M2; M4-M5..... 79-81  
 Midmark Synthesis Glove/ Towel/ Cup Dispenser D1-D4..... 82  
 Midmark Synthesis Dressing Nook N1-N2 ..... 83  
 Midmark Synthesis Waste Receptacle W1 ..... 84  
 Midmark Synthesis Instrument Panel I1-I2 ..... 85  
 Ritter 261 Waste Receptacle..... 86

**Instrument Processing**

Ritter M11 Steam Sterilizer ..... 87  
 Ritter M9 Steam Sterilizer ..... 87  
 Ritter M9D Manual Sterilizer..... 88  
 Ritter M3 Steam Sterilizer ..... 88  
 Sterilizer Accessories..... 89  
 VistaCool™ Direct-to-Drain System ..... 90  
 QuickClean™ Ultrasonic Cleaner QC6/QC6R..... 91  
 QuickClean™ Ultrasonic Cleaner QC3/QC3R..... 91  
 QuickClean™ Ultrasonic Cleaner QC1 ..... 92  
 Ultrasonic Parts and Accessories..... 92

## Midmark Customer Experience

To place an order, email: [Orders-Medical@midmark.com](mailto:Orders-Medical@midmark.com). For order status, requests and information, contact Medical Customer Experience at: **1.844.856.1233, Option 1**  
**Email: [allcustservice-medical@midmark.com](mailto:allcustservice-medical@midmark.com)**

## Midmark Technical Service

For service part identification, technical/operational questions, warranty and product issues/complaints, contact Technical Service at: **1.844.856.1230**.

**Sterilizer + Ultrasonic Products:** [sterilizers.ts@midmark.com](mailto:sterilizers.ts@midmark.com)

**Diagnostic Products:** [techsupport@midmark.com](mailto:techsupport@midmark.com)

**Procedure + Exam Chairs, Seating, Lighting, Diagnostic Instrument Wall/Desk Kits, Otoscopies,**

**Ophthalmoscopes + Sphygmomanometers:** [examroomeq.ts@midmark.com](mailto:examroomeq.ts@midmark.com)

**Medical Cabinetry, Carts + Workstations:** [casework.ts@midmark.com](mailto:casework.ts@midmark.com)

## Special Notes

1. **ORDERING**—Please provide the following information on all orders:

- Purchase order number
- Bill to and ship to addresses
- Product model and/or part number(s)
- Description including upholstery color and/or panel/countertop color
- Quantity needed
- Pricing
- Required shipping date
- Special shipping requirements (see details below)  
MINIMUM ORDER — \$100.

### DROP SHIP REQUIREMENTS

- Please specify if the customer has a loading dock or has special shipping requirements.
- Receiving contact name, email and phone number is required.

2. **SPECIAL DESIGN** — Orders for equipment of special design, construction, upholstery or special finish are not subject to cancellation after your acceptance of our acknowledgment.

3. **SPECIAL SHIPPING REQUIREMENTS** — These services are available to enhance delivery to your customers. Please call for pricing.

- **CALL BEFORE:** Carrier will call delivery contact to provide a drop-off timeframe once they have possession of the product from Midmark and are at the final terminal.
- **LIFTGATE:** Carrier uses liftgate/ramp to get product to ground level when facility does not have a dock. Facility is responsible for taking product once at ground level.
- **INSIDE DELIVERY:** Carrier unloads product and delivers inside the first set of doors only, provided entrance to facility can accommodate freight as packaged (including skids) and there are no steps.
- **MIDMARK DELIVERY SERVICE:** Carrier sets product in place per the delivery contact. Carrier is responsible for assembling product and removal of all packaging, trash and skids from the delivery site (custom casework does not include installation).

4. **FREIGHT DAMAGE** — Merchandise arriving with external damage should be so described on the bill of lading. In case of concealed damage, carrier must be notified within 15 days after receipt and claims for damage instituted.

5. **ORDER CHANGES** — Any changes made to an order within 48 hours of the scheduled ship date may affect the scheduled ship date, requiring Midmark to reschedule.

6. **ORDER CANCELLATION** — Any order, if cancelled within ten days before the scheduled shipping date, will be charged a 20% cancellation charge.

7. **RETURN OF MERCHANDISE** — Authorization must be obtained for the return of any item. If return is authorized, transportation charges must be prepaid. Risk of loss with respect to returned goods shall remain with purchaser until receipt by Midmark. Items returned for credit are subject to a charge of 25% or more to cover inspection, damage in transit, reconditioning and repacking. Items may be returned within 6 months (except service parts within 90 days) from the original ship date if in a good, unused condition. RMA # should be written on the outside of the packaging. Diagnostic products may be returned within 30 days if in a good, unused condition and in original sealed packaging.

**Note: Certain items are non-returnable: Specials, Synthesis® Cabinetry, Midmark Workstations, Discontinued Items or Parts, Custom or Non-standard Colors, Service Upholstery and Service Parts if value is less than \$50. Diagnostic Consumables e.g. Electrodes, Prep Kits, Pneumotach Mouthpieces, Temperature Probe Covers, Nose Clips, Thermal Paper, All Disposable Items and Diagnostic Items value less than \$50 per.**

8. **SPECIAL ORDER REQUIREMENTS** — Any order requiring any special deviation from standard products will carry the following disclaimer:

We hereby accept and acknowledge the above captioned order, provided however, because of the special nature of this order in that this equipment is being supplied upon special design specifications, special materials and/or field modifications, as per your request, the standard warranties of Midmark and Ritter on its products are not applicable to this order, and this acceptance and acknowledgment is conditioned upon your indemnification and promise to hold Midmark and Ritter, its officers, employees and agents, harmless in respect to any and all warranty claims, damages, both direct and consequential, expenses, costs and any liability arising from the use of or relating to this special equipment.

9. **IMPROVE PRODUCT QUALITY** — Midmark Corporation reserves the right to change design and construction without notice.

## Midmark and Ritter Medical Equipment (Tables, Chairs, Physician Seating, Lights, Workstations and Sterilizers) Limited Warranty

**SCOPE OF WARRANTY** Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace components of the domestic and international medical products manufactured by Midmark (except for components not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

Midmark warrants to the original retail purchaser that during the applicable warranty period it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications.

**APPLICABLE WARRANTY PERIOD** The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product and shall be one (1) year for all warranted products and components.

**OBTAINING WARRANTY SERVICE** Warranty service must be obtained through either Midmark or an authorized dealer in the Midmark product line for which warranty service is requested. Midmark may be contacted for warranty service inquiries or issues via email at [www.midmark.com](mailto:www.midmark.com); by phone at 1-800-MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380.

It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

**EXCLUSIONS:** This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Products considered to be of a consumable or sterile nature;
- (4) Accessories or parts not manufactured by Midmark;
- (5) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (6) Costs and expenses of routine maintenance and cleaning; and
- (7) Representations and warranties made by any person or entity other than Midmark.
- (8) Matching of color, grain or texture except to commercially acceptable standards;
- (9) Changes in color caused by natural or artificial light;
- (10) Custom manufactured products;
- (11) Alterations or modifications to the product by any person or entity other than Midmark; and
- (12) Products that would otherwise be covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market.

**SOFTWARE; WITH RESPECT TO SOFTWARE THAT IS A PRODUCT OR COMPONENT THEREOF, MIDMARK DOES NOT WARRANT THAT THE SOFTWARE: (1) IS ERROR FREE; (2) CAN BE USED WITHOUT PROBLEMS OR INTERRUPTIONS; OR (3) IS FREE FROM VULNERABILITY TO INTRUSION OR ATTACK BY VIRUSES OR OTHER METHODS.**

**EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER** MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN.

**WARRANTY DISCLAIMER** THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

**STATUTE OF LIMITATIONS** No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

**NO AUTHORIZATION** No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

Warranty information included in product manuals



# Midmark 641

## Procedure Chair

- Power height, back, foot and tilt
- 450-pound patient weight capacity
- 19" to 40" seat height adjustment (non-rotation models)
- 22.5" to 43.5" seat height adjustment (rotation models)
- Perimeter contouring for maximum access to patient
- Flip-out foot extension
- Hand and foot control standard
- Paper roll holder standard
- Choice of Premium/Premium UltraFree\* seamless upholstery
- Meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards
- Table weight: 450 lb
- Shipping weight: 500 lb

### Model

### Description

#### Standard Models

641-002 Procedure Chair Non-Programmable w/o receptacle, w/o top  
002-0966-00-XXX Upholstery top - Premium

641-002 Procedure Chair Non-Programmable w/o receptacle, w/o top  
002-10147-XXX 641 Premium UltraFree\* 24" Upholstery

641-003 Procedure Chair Programmable w/receptacle, w/o top  
002-0966-00-XXX Upholstery top - Premium

641-003 Procedure Chair Programmable w/o receptacle, w/o top  
002-10147-XXX 641 Premium UltraFree\* 24" Upholstery

#### Rotation Models

641-004 Procedure Chair Non-Programmable w/rotation, w/o top  
002-0966-00-XXX Upholstery top - Premium

641-004 Procedure Chair Non-Programmable w/rotation, w/o top  
002-10147-XXX 641 Premium UltraFree\* 24" Upholstery

641-005 Procedure Chair Programmable w/rotation, w/o top  
002-0966-00-XXX Upholstery top - Premium

641-005 Procedure Chair Programmable w/rotation, w/o top  
002-10147-00-XXX 641 Premium UltraFree\* 24" Upholstery

Accessories	Description
<b>Arms/Armboards</b>	
9A179002-XX <sup>7</sup>	Fixed armboard
9A396001-XXX <sup>5,7</sup>	Chair arms
9A396001-XXX	Description: Chair arms, UltraFree*
9A81001-XXX <sup>7</sup>	Articulating armboard
9A82001-XXX <sup>7</sup>	Hand surgery armboard
<b>Casters</b>	
9A330001 <sup>3</sup>	Caster base, swivel wheels
9A365001 <sup>3</sup>	Casters
<b>Controls</b>	
9A408001 <sup>4</sup>	Wireless control, non-programmable, linear-style, plant installed
9A408002 <sup>4</sup>	Wireless control, programmable, linear-style, plant installed
9A408003 <sup>4</sup>	Wireless control, non-programmable, linear-style, field installed
9A408004 <sup>4</sup>	Wireless control, programmable, linear-style, field installed
<b>Headrests</b>	
9A378001-XXX <sup>7</sup>	Flat headrest, Premium upholstery
9A378001-XXX	Flat headrest, UltraFree* upholstery
9A382001-XXX <sup>7</sup>	U-shaped headrest, Premium upholstery
9A382001-XXX	U-shaped headrest, UltraFree* upholstery
9A386001-XXX <sup>7</sup>	Oval-shaped headrest, Premium upholstery
9A386001-XXX	Oval-shaped headrest, UltraFree* upholstery
9A390001-XXX <sup>7</sup>	Rectangular-shaped headrest, Premium upholstery
9A390001-XXX	Rectangular-shaped headrest, UltraFree* upholstery
9A394001-XXX <sup>7</sup>	Magnetic headrest, Premium upholstery
9A394001-XXX	Magnetic headrest, UltraFree* upholstery
<b>Instrument Trays/Work Shelves/IV Poles</b>	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A267003 <sup>6</sup>	Instrument tray, double arm
9A77001 <sup>6</sup>	IV pole
<b>Miscellaneous</b>	
9A152010	3' power cord
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A406001	Security straps, plant installed
9A406002	Security straps, field installed
9A78001-XXX <sup>6,7</sup>	Vision block screen
9A416001	Foot section plastic cover
002-0966-00-999	Upholstery top - Premium Special Color
002-0967-00-XXX	CAL 133 Upholstery top - Premium
002-0967-00-999	CAL 133 Upholstery top - Premium Special Color
002-10148-XXX	CAL 133 Upholstery Top - Premium UltraFree*

<sup>3</sup> - Not available on rotation models

<sup>4</sup> - Accessory supplied with hand and foot control (batteries included)

<sup>5</sup> - Interference with armboard accessories will likely occur

<sup>6</sup> - Requires 9A20500X base rail

<sup>7</sup> - Additional charge for Special, Premium UltraFree\* for Accessory

\*ULTRALEATHER® by Ultrafabrics®



# Midmark 630

## HUMANFORM® Procedure Chair

- Power height, back, foot and tilt
- Motion Profiling® designed to provide more natural, human-like movement that works with the patient and provider for a clear view of the procedure site and improved comfort
- Active Sensing Technology® (AST) is designed to actively protect patient and provider by sensing and avoiding harmful impact during movement, such as objects being caught under the seat and footrest
- Barrier-Free design
- 650-pound patient weight capacity
- 18" to 40" seat height adjustment (non-rotational models)
- 21.5" to 43.5" seat height adjustment (rotational models)
- Hand and foot controls standard
- Available in 28" or 32" upholstery widths
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Table weight: 450 lb
- Shipping weight: 525 lb

Model	Description
<b>Non-Programmable</b>	
630-010 002-10108-XXX	Procedure Chair Non-Programmable w/receptacle, w/o top 630 Premium 28" Upholstery
630-010 002-10114-XXX	Procedure Chair Non-Programmable w/receptacle, w/o top 630 Premium UltraFree* 28" Upholstery
630-011 002-10108-XXX	Procedure Chair Non-Programmable wireless, w/receptacle, w/o top 630 Premium 28" Upholstery
630-011 002-10114-XXX	Procedure Chair Non-Programmable wireless, w/receptacle, w/o top 630 Premium UltraFree* 28" Upholstery
<b>Non-Programmable with Rotation</b>	
630-012 002-10108-XXX	Procedure Chair Non-Programmable w/rotation, w/o top 630 Premium 28" Upholstery
630-012 002-10114-XXX	Procedure Chair Non-Programmable w/rotation, w/o top 630 Premium UltraFree* 28" Upholstery
630-013 002-10108-XXX	Procedure Chair Non-Programmable wireless, w/rotation, w/o top 630 Premium 28" Upholstery
630-013 002-10114-XXX	Procedure Chair Non-Programmable wireless, w/rotation, w/o top 630 Premium UltraFree* 28" Upholstery



Model	Description
<b>Programmable</b> 630-020 002-10108-XXX	Procedure Chair Programmable w/receptacle, w/o top 630 Premium 28" Upholstery
630-020 002-10114-XXX	Procedure Chair Programmable, w/receptacle, w/o top 630 Premium UltraFree* 28" Upholstery
630-021 002-10108-XXX	Procedure Chair Programmable wireless, w/receptacle, w/o top 630 Premium 28" Upholstery
630-021 002-10114-XXX	Procedure Chair Programmable wireless, w/receptacle, w/o top 630 Premium UltraFree* 28" Upholstery
<b>Programmable with Rotation</b>	
630-022 002-10108-XXX	Procedure Chair Programmable, w/rotation, w/o top 630 Premium 28" Upholstery
630-022 002-10114-XXX	Procedure Chair Programmable, w/rotation, w/o top 630 Premium UltraFree* 28" Upholstery
630-023 002-10108-XXX	Procedure Chair Programmable wireless, w/rotation, w/o top 630 Premium 28" Upholstery
630-023 002-10114-XXX	Procedure Chair Programmable wireless, w/rotation, w/o top 630 Premium UltraFree* 28" Upholstery
<b>Upholstery Options</b>	
002-10108-999 002-10114-999	Upholstery top - 28" wide Premium Special Color Upholstery top - 28" wide UltraFree* Special Color
002-10112-XXX 002-10112-999	Upholstery top - 32" wide Premium Upholstery top - 32" wide Premium Special Color
002-10116-XXX 002-10116-999	Upholstery top - 32" wide UltraFree* Upholstery top - 32" wide UltraFree* Special Color
002-10109-XXX 002-10109-999	CAL 133 Upholstery top - 28" wide Premium CAL 133 Upholstery top - 28" wide Premium Special Color
002-10115-XXX 002-10115-999	CAL 133 Upholstery top - 28" wide UltraFree* CAL 133 Upholstery top - 28" wide UltraFree* Special Color
002-10113-XXX 002-10113-999	CAL 133 Upholstery top - 32" wide Premium CAL 133 Upholstery top - 32" wide Premium Special Color
002-10117-XXX 002-10117-999	CAL 133 Upholstery top - 32" wide UltraFree* CAL 133 Upholstery top - 32" wide UltraFree* Special Color

\*ULTRALEATHER® by Ultrafabrics®

Accessories	Description
<b>Arm Systems/Armboards</b>	
9A550001-XXX <sup>7</sup>	Chair arms, 28" wide Premium upholstery
9A550001-XXX	Chair arms, 28" wide UltraFree* upholstery
9A550005-XXX <sup>7</sup>	Chair arms, 32" wide Premium upholstery
9A550005-XXX	Chair arms, 32" wide UltraFree* upholstery
9A179002-XXX <sup>6,7</sup>	Fixed Armboard
9A81001-XXX <sup>6,7</sup>	Articulating Armboard
9A82001-XXX <sup>6,7</sup>	Hand surgery Armboard
<b>Controls (additional)</b>	
9A552001	Wired hand control, non-programmable, field installed
9A553001	Wired foot control, non-programmable, field installed
9A552002	Wired hand control, programmable, field installed
9A553002	Wired foot control, programmable, field installed
<b>Headrests</b>	
9A496001-XXX <sup>7</sup>	Rectangular-shaped headrest, Premium upholstery
9A496001-XXX	Rectangular-shaped headrest, UltraFree* upholstery
9A500001-XXX <sup>7</sup>	Oval-shaped headrest, Premium upholstery
9A500001-XXX	Oval-shaped headrest, UltraFree* upholstery
9A504001-XXX <sup>7</sup>	Round headrest, Premium upholstery
9A504001-XXX	Round headrest, UltraFree* upholstery
9A508001-XXX <sup>7</sup>	Flat headrest, Premium upholstery
9A508001-XXX	Flat headrest, UltraFree* upholstery
9A512001-XXX <sup>7</sup>	U-shaped headrest, Premium upholstery
9A512001-XXX	U-shaped headrest, UltraFree* upholstery
9A516001-XXX <sup>7</sup>	Magnetic headrest, Premium upholstery
9A516001-XXX	Magnetic headrest, UltraFree* upholstery
<b>Instrument Trays/Work Shelves/IV Poles</b>	
9A267003 <sup>4</sup>	Instrument tray, double arm
9A77001 <sup>4</sup>	IV pole
<b>Knee Crutches</b>	
9A410001 <sup>5</sup>	Standard knee crutches (crutches only)
9A410002 <sup>5</sup>	Standard knee crutch system
9A411001 <sup>5</sup>	Articulating knee crutches (crutches only)
9A411002 <sup>5</sup>	Articulating knee crutch system
9A412001 <sup>5</sup>	Knee crutch mounting brackets (brackets only)
<b>Accessory Rails</b>	
9A521001	Seat Accessory Rails, Surgical Size, 28" wide upholstery, plant installed
9A522001	Seat Accessory Rails, Surgical Size, 32" wide upholstery, plant installed
9A567001	Seat Accessory Rails, Surgical Size, 28" wide upholstery, field installed
9A568001	Seat Accessory Rails, Surgical Size, 32" wide upholstery, field installed
9A557001	Accessory Rails, 28" wide upholstery, field installed
9A558001	Accessory Rails, Surgical Size, 28" wide upholstery, field installed
9A557002	Accessory Rails, 32" wide upholstery, field installed
9A558002	Accessory Rails, Surgical Size, 32" wide upholstery, field installed
<b>Miscellaneous</b>	
9A554001 <sup>3</sup>	Casters - set of 4
9A152010	3' power cord
9A157003 <sup>4</sup>	Security side panels
9A559001 <sup>6</sup>	Security straps, field installed
9A555001	Gynecology debris tray accessory
9A556001	Gynecology debris tray (tray only)
9A78001-XXX <sup>7</sup>	Anesthesia/Vision block screen
AD 1779 <sup>3</sup>	Tamper-resistant duplex receptacle

<sup>3</sup> - Not available on rotation models

<sup>4</sup> - Requires (1) 9A557001/9A557002/9A558001/9A558002 accessory rails

<sup>5</sup> - Not compatible with 9A521001, 9A522001, 9A567001 or 9A568001

<sup>6</sup> - Requires (1) 9A557001/9A557002/9A558001/9A558002 accessory rails and (1) 9A521001/9A522001 seat rails

<sup>7</sup> - Additional charge for Special, UltraFree\* Color for Accessory

\*ULTRALEATHER® by Ultrafabrics®



# Midmark 647

## Podiatry Chair

- Power height, tilt and back
- 450-pound patient weight capacity
- 19" to 40" seat height adjustment (non-rotation models)
- 22.5" to 43.5" seat height adjustment (rotation models)
- 9" foot section extension
- Sealed foot control standard
- Integrated touch controls standard
- Chair arms standard
- Choice of Premium seamless UltraFree upholstery
- Protective foot cover standard with upholstery
- Meets UL 60601-1 and CAN/CSA C22.2 No. 601.1-M90 standards
- Table weight: 433 lb (non-rotation models)
- Shipping weight: 515 lb (non-rotation models)
- Table weight: 539 lb (rotation models)
- Shipping weight: 621 lb (rotation models)

Model	Description
647-001 002-1119-XXX	Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Premium
647-001 002-10129-XXX	Podiatry Chair Non-programmable w/receptacles, w/o top 646/647 Premium UltraFree* 28" Upholstery
647-002 002-1119-XXX	Podiatry Chair Programmable w/receptacles, w/o top Upholstery top - Premium
647-002 002-10129-XXX	Podiatry Chair Programmable w/receptacles, w/o top 646/647 Premium UltraFree* 28" Upholstery
647-003 002-1119-XXX	Podiatry Chair Non-programmable w/rotation, w/o receptacles & top Upholstery top - Premium
647-003 002-10129-XXX	Podiatry Chair Non-programmable w/rotation, w/o receptacles & top 646/647 Premium UltraFree* 28" Upholstery
647-004 002-1119-XXX	Podiatry Chair Programmable w/rotation, w/o receptacles & top Upholstery top - Premium
647-004 002-10129-XXX	Podiatry Chair Programmable w/rotation, w/o receptacles & top 646/647 Premium UltraFree* 28" Upholstery

\*ULTRALEATHER® by Ultrafabrics®

Accessories	Description
<b>Instrument Trays/Work Shelves/IV Poles</b>	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A432001	Instrument tray, swing arm (tray only)
9A77001 <sup>4</sup>	IV Pole
<b>Casters</b>	
9A330001 <sup>3</sup>	Caster base, swivel wheels
9A365001 <sup>3</sup>	Casters
<b>Miscellaneous</b>	
9A418001	Wireless foot control, non-programmable, linear-style, plant installed
9A418002	Wireless foot control, non-programmable, linear-style, field installed
9A418003	Wireless foot control, programmable, linear-style, plant installed
9A418004	Wireless foot control, programmable, linear-style, field installed
9A422001-XXX <sup>5</sup>	Vision block screen
9A420001	Plastic protective foot cover
9A423001	Paper roll holder
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A152010	3' power cord
002-1119-999	Upholstery top - Premium Special Color
002-1120-XXX	CAL 133 Upholstery top - Premium
002-1120-999	CAL 133 Upholstery top - Premium Special Color
002-10130-XXX	CAL 133 Upholstery top - UltraFree*

<sup>3</sup> - Not available on rotation models

<sup>4</sup> - Requires 9A20500X base rail

<sup>5</sup> - Additional charge for Special or UltraFree\* Accessory

\*ULTRALEATHER® by Ultrafabrics®



# Midmark 646

## Podiatry Chair

Not available for  
sale in Canada

- Power tilt and back
- 450-pound patient weight capacity
- 21" fixed seat height
- 9" foot section extension
- Sealed foot control standard
- Integrated touch controls standard
- Chair arms standard
- Choice of Premium Seamless or UltraFree\* upholstery
- Protective foot cover standard with upholstery
- Meets UL 60601-1
- Table weight: 323 lb
- Shipping weight: 405 lb

Model	Description
646-001 002-1119-XXX	Podiatry Chair Non-programmable w/receptacles, w/o top Upholstery top - Premium
646-001 002-10129-XXX	Podiatry Chair Non-programmable w/ receptacles, w/o top 646/647 Premium UltraFree* 28in Upholstery

Accessories	Description
<b>Instrument Trays/Work Shelves/IV Poles</b>	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A432001	Instrument tray, swing arm (tray only)
9A77001 <sup>3</sup>	IV Pole
<b>Casters</b>	
9A330001	Caster base, swivel wheels
9A365001	Casters
<b>Miscellaneous</b>	
9A422001-XXX <sup>4</sup>	Vision block screen
9A420001	Plastic protective foot cover
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A152010	3' power cord
002-1119-999	Upholstery top - Premium Special Color
002-1120-XXX	CAL 133 Upholstery top - Premium
002-1120-999	CAL 133 Upholstery top - Premium Special Color

<sup>3</sup> - Requires 9A20500X base rail

<sup>4</sup> - Additional charge for Special or UltraFree\* Upholstery Color for Accessory

\*ULTRALEATHER® by Ultrafabrics®



## Midmark 626 Barrier-Free® Examination Chair

- 650-pound weight capacity
- Low height of 15.5" (uncompressed) to 37" high
- Clean Assist Roller System
- Power base and back
- Hand and foot controls
- Digital Scale (-003 and -004 models)
- Heated upholstery options (-002 and -004 models)
- Quick Exam® and Return to Home
- Active Sensing Technology®
- Mist and Stone table base color
- Lockout of hand and foot controls
- Table weight: 435 lb
- Shipping weight: 485 lb

Model	Description
<b>Standard Models</b>	
626-001 002-2000-XXX <sup>1</sup>	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, w/o top 626 Premium 28" Upholstery
626-001 002-2002-XXX <sup>1</sup>	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, w/o top 626 Premium 32" Upholstery
626-001 002-2004-XXX <sup>1</sup>	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, w/o top 626 Premium UltraFree* 28" Upholstery
626-002 002-2006-XXX <sup>1</sup>	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, heated upholstery options, w/o top 626 Premium 28" Heated Upholstery
626-002 002-10201-XXX <sup>1</sup>	Barrier-Free Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, heated upholstery options, w/o top 626 Premium 32" Heated Upholstery

\*ULTRALEATHER® by Ultrafabrics®

Model	Description
<b>Standard Models</b>	
626-002	Barrier-Free® Examination Chair w/ Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, wired controls, heated upholstery options, w/o top
002-2011-XXX <sup>1</sup>	626 Premium UltraFree Heated 28" Upholstery
626-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2000-XXX <sup>1</sup>	626 Premium 28" Upholstery
626-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2002-XXX <sup>1</sup>	626 Premium 32" Upholstery
626-003	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2004-XXX <sup>1</sup>	626 Premium UltraFree* 28" Upholstery
626-004	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2006-XXX <sup>1</sup>	626 Premium 28" Heated Upholstery
626-004	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-10201-XXX <sup>1</sup>	626 Premium 32" Heated Upholstery
626-004	Barrier-Free Examination Chair w/ Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wireless controls, w/o top
002-2011-XXX <sup>1</sup>	626 Premium UltraFree* Heated 28" Upholstery

\*ULTRALEATHER® by Ultrafabrics®

<b>Model</b>	<b>Description</b>
<b>Standard Models</b>	
626-005	Barrier-Free Examination Chair w/Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2000-XXX <sup>1</sup>	626 Premium 28" Upholstery
626-005	Barrier-Free Examination Chair w/Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2002-XXX <sup>1</sup>	626 Premium 32" Upholstery
626-005	Barrier-Free Examination Chair w/Clean Assist Roller System, Digital Scale, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2004-XXX <sup>1</sup>	626 Premium UltraFree* 28" Upholstery
626-006	Barrier-Free Examination Chair w/Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2006-XXX <sup>1</sup>	626 Premium 28" Heated Upholstery
626-006	Barrier-Free Examination Chair w/Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-10201-XXX <sup>1</sup>	626 Premium 32" Heated Upholstery
626-006	Barrier-Free Examination Chair w/Clean Assist Roller System, Digital Scale, heated upholstery options, receptacles, drawer heater, pelvic tilt, wired controls, w/o top
002-2011-XXX <sup>1</sup>	626 Premium UltraFree* 28" Heated Upholstery



## Arm Systems

Part #	Description
9A600002 <sup>2</sup>	Patient Support Rails Plus, field installed
9A600001 <sup>2</sup>	Patient Support Rails, field installed
9A582001	Patient Support Rails, receiver extension - pair; field install for 32" wide upholstery

## Knee Crutches

9A643001	Standard knee crutch kit w/ acc. receiver, field install
9A645001	Articulating knee crutch kit w/ acc. receiver, field install
9A641001	Mounting brackets, hardware w/ acc. receiver, field install

## Vital Signs Accessories

9A478001	IQscale serial cable for IQvitals Vital Signs, 6 foot length, coiled
9A478002	IQscale serial cable for IQvitals Vital Signs, 15 foot length, straight
9A478003	IQscale serial cable for IQvitals Vital Signs, 30 foot length, straight
9A478004	IQscale serial cable for IQvitals Vital Signs, 50 foot length, straight
9A629001	IQvitals Vital Signs table mount with USB and serial cables
9A629002	IQvitals Vital Signs table mount with serial cable
9A629003	IQvitals Vital Signs table mount, non-hub
9A479001	Welch Allyn Spot Vitals Lxi, 6 foot length coiled
9A479002	Welch Allyn Spot Vitals Lxi, 15 foot length straight
002-10064-00	626 chair connectivity kit

## Miscellaneous

9A677001	Wireless hand and foot controls
9A610001	Stainless steel treatment pan
002-2015-XXX <sup>1</sup>	626 Premium UltraFree* 32" upholstery
002-10203-XXX <sup>1</sup>	626 Premium UltraFree* heated 32" upholstery
002-2003-XXX <sup>1</sup>	626 Premium 32" upholstery, CAL 133
002-2001-XXX <sup>1</sup>	626 Premium 28" upholstery, CAL133
002-2005-XXX <sup>1</sup>	626 Premium UltraFree* 28" upholstery, CAL133
002-2016-XXX <sup>1</sup>	626 Premium UltraFree* 32" upholstery, CAL133
002-2007-XXX <sup>1</sup>	626 Premium 28" heated upholstery, CAL133
002-10202-XXX <sup>1</sup>	626 Premium 32" heated upholstery, CAL133
002-2012-XXX <sup>1</sup>	626 Premium UltraFree* heated 28" upholstery, CAL133
002-10204-XXX <sup>1</sup>	626 Premium UltraFree* heated 32" upholstery, CAL133
002-10360-XXX	224/225/626 Flat 28" upholstery, CAL133
002-10362-XXX	224/225/626 Flat UltraFree* 28" upholstery, CAL133

<sup>1</sup> - See color palette on page 19

<sup>2</sup> - 9A582001 receiver extensions required for use with 32" wide upholstery

\*ULTRALEATHER® by Ultrafabrics®




## Midmark 640 Pediatric Examination Table

- Pebble grey base and panel color only
- An upholstered end panel and beveled sides designed to help keep infants from rolling off
- Flat sides with no sharp edges
- Accurate and easy-to-use, built-in scale
- 36-inch exam surface height
- Ample storage and a hidden paper roll holder
- Meets ANSI/AAMI ES 60601-1:2005 & CAN/CSA C22.2 No. 60601-1-08 standards
- 42" W x 51.5" H x 23"D
- Table weight: 170 lb
- Shipping weight: 210 lb

Model	Description
640-001-XXX	Pediatric Examination Table w/digital scale
640-001-999	Pediatric Examination Table w/digital scale - special color
<b>Miscellaneous</b>	
640-011-XXX	CAL 133 Pediatric Examination Table w/digital scale
640-012-XXX	CAL 133 Pediatric Examination Table with door and drawer lock, plant installed
640-011-999	CAL 133 Pediatric Examination Table w/digital scale - special color
640-012-999	CAL 133 Pediatric Examination Table w/door and drawer lock, plant installed - special color
640-002-XXX	Pediatric Examination Table w/door and drawer lock, plant installed
640-002-999	Pediatric Examination Table w/door and drawer lock, plant installed - special color

# Upholstery Color Palette

XXX - XXX - XXX



Example:  
630-010 *Chair*  
002-10108-843 *Top*  
630 HUMANFORM® *Procedure Chair - Oasis, Premium Upholstery Top*

## Standard Upholstery Colors

Mist-813	Soothing Blue-855	Curative Copper-846
Stone-814	Dream-847	Cranberry-859
Restful Path-858	Mineral-844	Dark Linen-856
Iron Ore-851	Healing Waters-848	Latte-860
Lunar Gray-845	Shaded Garden-853	Robust Brown-852
Obsidian-857	Citrus-849	Deep Earth-854
Oasis-843	Sandy Retreat-850	

## UltraFree\* Upholstery Colors

UltraFree Stone-865	UltraFree Spice-863
UltraFree Restful Path-868	UltraFree Cranberry-869
UltraFree Obsidian-867	UltraFree Dark Linen-866
UltraFree Harbor-861	UltraFree Latte-870
UltraFree Wheatgrass-862	UltraFree Branch-864

\*ULTRALEATHER® by Ultrafabrics®



# Ritter 230

## Procedure Chair

- Power height, back, foot and tilt
- 450-pound patient weight capacity
- 22.5" to 40" seat height adjustment (non-rotation models)
- 26" to 43.5" seat height adjustment (rotation models)
- Knee-chest configurable
- Hand control standard
- Meets UL60601-1 and CAN/CSA C22.2 No. 601.1 – M90 standards
- Table weight: 455 lb
- Shipping weight: 525 lb

Model	Description
230-001 002-0824-XXX	Procedure Chair w/o top 230 Seamless Upholstery top
230-001 002-10133-XXX	Procedure Chair w/o top 230 Seamless UltraFree* 28" Upholstery
230-002 002-0824-XXX	Procedure Chair w/receptacle, w/o top 230 Seamless Upholstery top
230-002 002-10133-XXX	Procedure Chair w/receptacle, w/o top 230 Seamless UltraFree* 28" Upholstery
230-003 002-0824-XXX	Procedure Chair w/rotation, w/o top 230 Seamless Upholstery top
230-003 002-10133-XXX	Procedure Chair w/rotation, w/o top 230 Seamless UltraFree* 28" Upholstery

\*ULTRALEATHER® by Ultrafabrics®

Accessories	Description
<b>Arm Systems/Armboards</b>	
9A179002-XXX <sup>6</sup>	Fixed armboard
9A328001	Adjustable arm system, plant installed
9A328002	Adjustable arm system, field installed
9A350001	Adjustable arm system (receiver only), plant installed
9A350002	Adjustable arm system (receiver only), field installed
9A350003	Adjustable arm system (arms only)
9A81001-XXX <sup>6</sup>	Articulating armboard
9A82001-XXX <sup>6</sup>	Hand surgery armboard
<b>Casters</b>	
9A330001 <sup>4</sup>	Caster base, swivel wheels
9A365001 <sup>4</sup>	Casters
<b>Controls</b>	
9A371003	Foot control, linear style
<b>Headrests/Footrests</b>	
9A337001-XXX <sup>6</sup>	Larger 24" wide footrest
9A51001-XXX <sup>6</sup>	Facial pad
<b>Instrument Trays/Work Shelves/IV Poles</b>	
9A427001	Swing arm instrument tray, fixed top
9A427002	Swing arm instrument tray, extension top
9A428001	Swing arm PC work shelf
9A267003 <sup>5</sup>	Instrument tray, double arm
9A77001 <sup>5</sup>	IV pole
<b>Knee Crutches</b>	
9A410001	Standard knee crutches (crutches only)
9A410002	Standard knee crutch system
9A411001	Articulating knee crutches (crutches only)
9A411002	Articulating knee crutch system
9A412001	Knee crutch mounting brackets (brackets only)
<b>Miscellaneous</b>	
9A152010	3' Power Cord
9A157003	Security side panels
9A205003	Base rails, plant installed
9A205004	Base rails, field installed
9A214001	Security straps, plant installed
9A214002	Security straps, field installed
9A329003 <sup>4</sup>	Urology drain system, plant installed
9A329004 <sup>4</sup>	Urology drain system, field installed
9A78001-XXX <sup>5,6</sup>	Vision block screen
002-0824-999	Seamless Upholstery top Special Color
002-0825-XXX	CAL 133 Upholstery Top-Seamless
002-0825-999	CAL 133 Upholstery Top - Seamless Special Color
002-0951-XXX	CAL 133 Flat Upholstery set, 24" width
002-0951-999	CAL 133 Flat Upholstery set, 24" width Special Color
002-0950-XXX	Flat upholstery set, 24" width
002-0950-999	Flat Upholstery set, 24" width Special Color
002-10135-XXX	Flat UltraFree* Upholstery
002-10134-XXX	CAL 133 Seamless UltraFree Upholstery
002-10136-XXX	CAL 133 Flat UltraFree* Upholstery

<sup>3</sup> - Not compatible with 9A337001

<sup>4</sup> - Not available on rotation models

<sup>5</sup> - Requires 9A20500X base rail

<sup>6</sup> - Additional charge for Special, UltraFree\* Accessory

\*ULTRALEATHER® by Ultrafabrics®



# Ritter 244

## Barrier-Free® Examination Table

- Power height, back and tilt
- Reverse trendelenburg
- 850-pound patient weight capacity
- 18" to 34" seat height adjustment
- 32" wide seat
- Hand control standard
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Table weight: 340 lb
- Shipping weight: 400 lb

Model	Description
244-001 002-0861-XXX	Barrier-Free Examination Table w/o top Seamless 32" Upholstery

244-001 002-10174-XXX	Barrier-Free Examination Table w/o top Seamless UltraFree* 32" Upholstery
--------------------------	--

### Accessories

9A627001	Adjustable arm system, field installed
9A346001	Casters
9A347001	Foot control
9A152010	3' Power Cord

### Miscellaneous

002-0861-999	Seamless 32" Upholstery Top Special Color
002-1017-XXX	CAL 133 Seamless 32" Upholstery Top
002-1017-999	CAL 133 Seamless 32" Upholstery Top Special Color
002-10175-XXX	CAL 133 Seamless UltraFree* 32" Upholstery Top

\*ULTRALEATHER® by Ultrafabrics®



# Ritter 225

## Barrier-Free® Examination Chair

- 500-pound weight capacity
- Low height of 18" to 37" high
- Power base and back
- Wired foot control
- Active Sensing Technology®
- Mist table base color table weight: 409 lb
- Shipping weight: 465 lb

Model	Description
<b>Standard Models</b>	
225-002	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top
002-2009-XXX <sup>1</sup>	224/225 Seamless 28" Upholstery
225-002	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top
002-10360-XXX <sup>1</sup>	224/225/626 Flat 28" Upholstery, CAL133
225-002	Barrier-Free Examination Chair w/ receptacles, drawer heater, pelvic tilt, w/o top
002-2013-XXX <sup>1</sup>	224/225 Seamless UltraFree* 28" Upholstery
225-003	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, w/o top
002-2009-XXX <sup>1</sup>	224/225 Seamless 28" Upholstery
225-003	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, w/o top
002-10360-XXX <sup>1</sup>	224/225/626 Flat 28" Upholstery, CAL133
225-003	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt, w/o top
002-2013-XXX <sup>1</sup>	224/225 Seamless UltraFree* 28" Upholstery
<b>Arm Systems</b>	
9A600002	Patient Support Rails Plus, Field Installed
9A600001	Patient Support Rails, Field Installed
9A615001	Assist Arms
<b>Knee Crutches</b>	
9A643001	Standard Knee Crutch Kit w/ Acc. Receiver, Field Install
9A645001	Articulating Knee Crutch Kit w/ Acc. Receiver, Field Install
9A641001	Mounting Brackets, Hardware w/ Acc. Receiver, Field Install
<b>Miscellaneous</b>	
9A610001	Stainless steel treatment pan
9A638001	Supply Assistant
9A681001	225 Hand Control Accessory
002-2010-XXX <sup>1</sup>	224/225 Seamless 28" Upholstery, CAL133
002-2014-XXX <sup>1</sup>	224/225 Seamless UltraFree* 28" Upholstery, CAL133
002-10362-XXX	224/225/626 Flat UltraFree* 28" Upholstery, CAL133
002-10149-XXX <sup>1</sup>	Premium Thick Footrest Pad, CAL 117/133
002-10214-XXX <sup>1</sup>	Premium UltraFree* Thick Footrest Pad, CAL 117/133

<sup>1</sup> - See color palette on page 19  
\*ULTRALEATHER® by Ultrafabrics®



# Ritter 224

## Barrier-Free® Examination Chair

- 500-pound weight capacity
- Low height of 18" to 37" high
- Power base
- Wired foot control
- Active Sensing Technology® (-002 and -003 models)
- Table weight: 409 lb
- Shipping weight: 465 lb

Model	Description
<b>Standard Models</b>	
224-001 002-2009-XXX <sup>1</sup>	Barrier-Free Examination Chair - excludes drawer and debris tray 224/225 Seamless 28" Upholstery
224-001 002-10360-XXX <sup>1</sup>	Barrier-Free Examination Chair - excludes drawer and debris tray 224/225/626 Flat 28" Upholstery, CAL133
224-001 002-2013-XXX <sup>1</sup>	Barrier-Free Examination Chair - excludes drawer and debris tray 224/225 Seamless UltraFree* 28" Upholstery
224-002 002-2009-XXX <sup>1</sup>	Barrier-Free Examination Chair w/receptacles, drawer heater, pelvic tilt, w/o top 224/225 Seamless 28" Upholstery
224-002 002-10360-XXX <sup>1</sup>	Barrier-Free Examination Chair w/receptacles, drawer heater, pelvic tilt, w/o top 224/225/626 Flat 28" Upholstery, CAL133
224-002 002-2013-XXX <sup>1</sup>	Barrier-Free Examination Chair w/receptacles, drawer heater, pelvic tilt, w/o top 224/225 Seamless UltraFree* 28" Upholstery
224-003 002-2009-XXX <sup>1</sup>	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt w/o top 224/225 Seamless 28" Upholstery
224-003 002-10360-XXX <sup>1</sup>	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt w/o top 224/225/626 Flat 28" Upholstery, CAL133
224-003 002-2013-XXX <sup>1</sup>	Barrier-Free Examination Chair w/Clean Assist Roller System, receptacles, drawer heater, pelvic tilt w/o top 224/225 Seamless UltraFree* 28" Upholstery

\*UltraFree by Ultrafabrics®



Model	Description
<b>Standard Models</b>	
<b>Arm Systems</b>	
9A600002	Patient Support Rails Plus, field installed
9A600001	Patient Support Rails, field installed
9A615001	Assist Arms
<b>Knee Crutches</b>	
9A643001	Standard knee crutch kit w/ acc. receiver, field install
9A645001	Articulating knee crutch kit w/ acc. receiver, field install
9A641001	Mounting brackets, hardware w/ acc. receiver, field install
<b>Miscellaneous</b>	
9A610001	Stainless steel treatment pan
9A638001	Supply Assistant
002-2010-XXX <sup>1</sup>	224/225 Seamless 28" Upholstery, CAL133
002-2014-XXX <sup>1</sup>	224/225 Seamless UltraFree* 28" Upholstery, CAL133
002-10362-XXX <sup>1</sup>	224/225/626 Flat UltraFree* 28" Upholstery, CAL133
002-10149-XXX <sup>1</sup>	Premium Thick Footrest Pad, CAL 117/133
002-10214-XXX <sup>1</sup>	Premium UltraFree* Thick Footrest Pad, CAL 117/133

<sup>1</sup> - See color palette on page 19



# Ritter 204

## Manual Examination Table

- 500-pound patient weight capacity
- Electrical receptacle standard
- Styled Seamless or Soft Touch stitched upholstery
- Pass-thru drawer design, no need to specify RH or LH
- Manual back with air-spring assist
- Meets UL 60601-1 and CAN/CSA C22.2 No.601.1-M90 standards
- Please specify upholstery/drawer panel color
- Mist table base color table weight: 250 lb
- Shipping weight: 300 lb

Model	Description
204-011 002-10137-XXX	Manual examination table w/o top 204 Soft Touch Uph, 28"
204-011 10138-XXX	Manual examination table w/o top 204 Upholstery top - Styled, Seamless
204-012 10137-XXX	Manual examination table w/drawer heater, pelvic tilt w/o top 204 Soft Touch Uph, 28"
204-012 10138-XXX	Manual examination table w/drawer heater, pelvic tilt w/o top 204 Upholstery top - Styled, Seamless

Accessories	Description
<b>Arm/Handle Systems</b>	
9A627001 <sup>3</sup>	Adjustable arm system (arms only)
9A632001 <sup>3</sup>	Assist handles (handles only), 28" wide upholstery only, field installed
9A632002	Accessory receiver only, 28" wide upholstery only, plant installed
9A632003	Accessory receiver only, 28" wide upholstery only, field installed
<b>Knee Crutches</b>	
9A410001	Standard knee crutch only (204/604/22X/62X/230/630) - field installed
9A410004	Standard knee crutch kit with accessory receiver (204/604/22X/62X) - field installed
9A411001	Articulating knee crutch only (204/604/22X/62X/230/630) - field installed
9A411004	Articulating knee crutch kit with accessory receiver (204/604/22X/62X) - field installed
9A412003	Mounting brackets, hardware with accessory receiver (204/604/22X/62X) - field installed
<b>Miscellaneous</b>	
9A104001	Urology drain pan w/hose
9A372001	Pass-through drawer stop (set of 3)
9A70000	Stainless steel treatment pan
9A629003	Midmark Digital Vital Signs table mount, non-hub
002-10143-XXX	204 Seamless UltraFree* Upholstery, 28"
002-10141-XXX	Flat Stitched Upholstery w/o pillow
002-10145-XXX	Flat UltraFree* Upholstery top w/o pillow
002-10193-00-XXX	Thick footrest pad kit (for Soft-Touch & Flat uph)
002-10194-00-XXX	Thick footrest pad kit (for Seamless/UltraFree* uph)

<sup>3</sup> - 9A632002 or 9A632003 accessory receiver required

\*ULTRALEATHER® by Ultrafabrics®



## Ritter 203

### Manual Examination Table

- 400-pound patient weight capacity
- 27" W x 74" D x 31" H
- Please specify upholstery/panel color
- Panels match upholstery color
- Table weight: 135 lb
- Shipping weight: 170 lb

Model	Description
203-011-XXX	Manual examination table standard
203-012-XXX	Manual examination table w/drawers
203-013-XXX	Manual examination table w/stirrups*
203-014-XXX	Manual examination table w/drawers and stirrups*
<b>Miscellaneous</b>	
203-011-999	Manual examination table standard - special color
203-012-999	Manual examination table w/drawers - special color
203-013-999	Manual examination table w/stirrups* - special color
203-014-999	Manual examination table drawers and stirrups* - special color

\* Stirrups factory installed only



## Ritter 95

### Manual Examination Table

- 350-pound patient weight capacity
- 27" W x 72" D x 30" H
- Please specify upholstery color
- Frame is a dark oak finish only
- Paper roll holder included
- Table weight: 80 lb
- Shipping weight: 100 lb

Model	Description
95-001-XXX <sup>3</sup>	Manual examination Table w/o shelf
95-001-999	Manual examination Table w/o shelf Special Color
<b>Accessories</b>	
9A364001-XXX <sup>3</sup>	<b>Description</b> Shelf (field installed only)
9A364001-999	Shelf Special Color
<b>Options</b>	
95-002-XXX	CAL 133 Manual Examination Table w/o shelf
95-002-999	CAL 133 Manual Examination Table w/o shelf - Special Color

<sup>3</sup> - Additional charge for UltraFree\* Upholstery/Accessory Color.



## Midmark 425

### Physician Seating (Hand Operated)

- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

#### Model

425-001  
002-0273-XXX<sup>1,2,3</sup>

#### Description

Basic stool w/o seat cushion  
Seat cushion

#### Accessories

002-0274-XXX<sup>1,2,4</sup>  
002-0275-XXX<sup>1,2,4</sup>  
002-1546-00

Back kit  
Arm kit  
Lab height w/black foot ring (425 only)

#### Miscellaneous

002-0944-XXX<sup>1,2</sup>  
002-0945-XXX<sup>1,2</sup>  
002-0946-XXX<sup>1,2</sup>

CAL 133 Seat Cushion-Premium Upholstery  
CAL 133 Back Kit-Premium Upholstery  
CAL 133 Arm Kit-Premium Upholstery

<sup>1</sup> - See color code chart on page 35

<sup>2</sup> - Additional charge for Special or UltraFree\* Upholstery Color

<sup>3</sup> - Additional charge for Special or UltraFree\* Upholstery Color for Accessory



## Midmark 427

### Physician Seating (Foot Operated)

- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

#### Model

427-001  
002-0273-XXX<sup>3</sup>

#### Description

Basic stool w/o seat cushion  
Seat cushion

#### Accessories

002-0274-XXX<sup>4</sup>  
002-0275-XXX<sup>4</sup>

Back kit  
Arm kit

#### Miscellaneous

002-0944-XXX  
002-0945-XXX  
002-0946-XXX

CAL 133 Seat Cushion-Premium Upholstery  
CAL 133 Back Kit-Premium Upholstery  
CAL 133 Arm Kit-Premium Upholstery

<sup>3</sup> - Additional charge for Special or UltraFree\* Upholstery Color

<sup>4</sup> - Additional charge for Special or UltraFree\* Upholstery Color for Accessory

\*UltraFree by Ultrafabrics®



## Ritter 279

### Physician Seating (Hand Operated)

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 30 lb
- Shipping weight: 35 lb

#### Model

279-001-XXX<sup>3</sup>

#### Description

Basic stool w/hand release - pneumatically adjustable



## Ritter 278

### Physician Seating (Foot Operated)

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 30 lb
- Shipping weight: 35 lb

#### Model

278-001-XXX<sup>3</sup>

#### Description

Basic stool w/foot release - pneumatically adjustable



hand release

## Ritter 277

### Physician Seating with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

#### Model

277-001-XXX<sup>3</sup>  
277-201-XXX<sup>3</sup>

#### Description

Basic stool w/back - pneumatically adjustable, hand release  
Basic stool w/back - pneumatically adjustable, foot release



foot release



hand release

## Ritter 276

### Physician Seating

- Chrome five leg base
- Weight capacity: 350 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

#### Model

276-001-XXX<sup>3</sup>  
276-201-XXX<sup>3</sup>

#### Description

Basic stool - pneumatically adjustable, hand release  
Basic stool - pneumatically adjustable, foot release



foot release

<sup>3</sup> - Additional charge for Special or UltraFree\* Upholstery Color



## Ritter 275

### Physician Seating with Backrest

- Chrome five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

**Model**  
275-001-XXX<sup>3</sup>

**Description**  
Basic stool w/back - manually adjustable



## Ritter 274

### Physician Seating

- Chrome five leg base
- Weight capacity: 350 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

**Model**  
274-001-XXX<sup>3</sup>

**Description**  
Basic stool - manually adjustable



## Ritter 273

### Physician Seating with Backrest

- Black composite five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

**Model**  
273-001-XXX<sup>3</sup>

**Description**  
Basic stool w/back - pneumatically adjustable



## Ritter 272

### Physician Seating

- Black composite five leg base
- Weight capacity: 300 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

**Model**  
272-001-XXX<sup>3</sup>

**Description**  
Basic stool - pneumatically adjustable

<sup>3</sup> - Additional charge for Special or UltraFree\* Upholstery Color

\*UltraFree by Ultrafabrics®



**Model**  
271-001-XXX<sup>3</sup>

## Ritter 271

### Physician Seating with Backrest

- Black composite five leg base
- Weight capacity: 225 lb
- Stool weight: 20 lb
- Shipping weight: 25 lb

#### Description

Basic stool w/back - manually adjustable



**Model**  
270-001-XXX<sup>3</sup>

## Ritter 270

### Physician Seating

- Black composite five leg base
- Weight capacity: 300 lb
- Stool weight: 15 lb
- Shipping weight: 20 lb

#### Description

Basic Stool - manually adjustable

#### Model

##### Miscellaneous

27X-X11-XXX

9A366001-XXX<sup>4</sup>

9A366002-XXX

27X-002-XXX<sup>3</sup>

27X-003-XXX<sup>3</sup>

27X-004-XXX<sup>3</sup>

27X-005-XXX<sup>3</sup>

9A367001

9A368001

#### Description

CAL 133 seat upholstery

Back rest kit (weight capacity, 225 lb)

CAL 133 back rest kit (weight capacity, 225 lb)

Set of soft rubber casters

Set of glides (not available on foot release models)

Set of locking casters

Auto locking casters

Foot ring, field installed (not available on foot release models)

6" lab height ext. kit (not available on foot release models)

<sup>3</sup> - Additional charge for Special or UltraFree\* Upholstery Color

<sup>4</sup> - Additional charge for Special or UltraFree\* Upholstery Color for Accessory

\*UltraFree by Ultrafabrics®





## Ritter 280

### Patient Side Chair

- 400-pound patient weight capacity
- Wall-saver, kickback rear leg design
- Leveling feet
- Chair weight: 25 lb
- Shipping weight: 35 lb

#### Model

280-001-XXX  
280-001-XXX  
280-002-XXX  
280-002-XXX

#### Description

Basic chair (without arms)  
Basic chair (without arms) UltraFree\* upholstery  
Basic chair (with arms)  
Basic chair (with arms) UltraFree\* upholstery

#### Miscellaneous

280-001-999 Basic chair (without arms) Special Color  
280-002-999 Basic chair (with arms) Special Color  
280-003-XXX CAL 133 Basic Chair (without arms)  
280-004-XXX CAL 133 Basic Chair (with arms)  
280-003-999 CAL 133 Basic Chair (without arms) Special Color  
280-004-999 CAL 133 Basic Chair (with arms) Special Color

See color code chart on page 35



## Ritter 282

### Bariatric Side Chair

- 650-pound patient weight capacity
- Wall-saver, kickback rear leg design
- Leveling feet
- Chair weight: 37 lb
- Shipping weight: 47 lb

#### Model

282-001-XXX  
282-001-XXX

#### Description

Basic chair  
Basic chair UltraFree\* upholstery

#### Miscellaneous

282-001-999 Basic chair Special Color  
282-002-XXX CAL 133 Basic Chair  
282-002-999 CAL 133 Basic Chair Special Color

See color code chart on page 35



## Midmark 680

Patient Side Chair (Discontinued 3/31/2019)

- 400-pound patient weight capacity
- Wall-saver, kickback rear leg design
- Leveling feet
- Chair weight: 20 lb
- Shipping weight: 30 lb

### Model

680-001-XXX  
680-001-XXX  
680-002-XXX  
680-002-XXX

### Description

Basic chair (without arms)  
Basic chair (without arms) UltraFree\* Upholstery  
Basic chair (with arms)  
Basic chair (with arms) UltraFree\* Upholstery

### Miscellaneous

680-001-999  
680-002-999  
680-003-XXX  
680-004-XXX  
680-003-999  
680-004-999

Basic chair (without arms) Special Color  
Basic chair (with arms) Special Color  
CAL 133 Basic Chair (without arms)  
CAL 133 Basic Chair (with arms)  
CAL 133 Basic Chair (without arms) Special Color  
CAL 133 Basic Chair (with arms) Special Color



## Ritter 281

Blood Draw Chair

- 400-pound patient weight capacity
- Optional storage drawer
- Chair weight: 60 lb
- Shipping weight: 75 lb



### Model

281-011-XXX  
281-011-XXX  
281-012-XXX  
281-012-XXX

### Description

Basic chair (without storage drawer)  
Basic chair (without storage drawer) UltraFree\* upholstery  
Basic chair (with storage drawer)  
Basic chair (with storage drawer) UltraFree\* upholstery


### Miscellaneous

281-011-999  
281-012-999  
281-013-XXX  
281-014-XXX  
281-013-999  
281-014-999

Basic chair (without storage drawer) Special Color  
Basic chair (with storage drawer) Special Color  
CAL 133 Basic Chair (without storage drawer)  
CAL 133 Basic Chair (with storage drawer)  
CAL 133 Basic Chair (without storage drawer) Special Color  
CAL 133 Basic Chair (with storage drawer) Special Color

\*UltraFree by Ultrafabrics®

# Seating Color-Specific Product Number Chart



Example: 276-001-814  
276 Physician Seating, Stone Upholstery

## Standard Upholstery Colors

Mist-813	Healing Waters-848
Stone-814	Shaded Garden-853
Restful Path-858	Citrus-849
Iron Ore-851	Sandy Retreat-850
Lunar Gray-845	Curative Copper-846
Obsidian-857	Cranberry-859
Oasis-843	Dark Linen-856
Soothing Blue-855	Latte-860
Dream-847	Robust Brown-852
Mineral-844	Deep Earth-854

## UltraFree\* Upholstery Colors

UltraFree Stone-865
UltraFree Restful Path-868
UltraFree Obsidian-867
UltraFree Harbor-861
UltraFree Wheatgrass-862
UltraFree Spice-863
UltraFree Cranberry-869
UltraFree Dark Linen-866
UltraFree Latte-870
UltraFree Branch-864

\*UltraFree by Ultrafabrics®

# Ritter 255

## LED Procedure Light

- 4,400 K color temperature
- Color Rendering Index (CRI) greater than 95
- 7,500 fc at 1 meter (39 inches) focal length
- Effortless positioning
- IEC/ES 60601-1, CAN/CSA C22.2 60601-1, 60601-2-41
- Weight: 65 lb (118 lb for dual light)
- Shipping weight: 120 lb (190 lb for dual light)



Mobile Procedure Light



Single Procedure Light



Dual Procedure Light



Wall Mount

Model	Description
255-001	Single mount, 8' ceiling
255-002	Single mount, 9' ceiling
255-003	Dual mount, 8' ceiling
255-004	Dual mount, 9' ceiling
255-005	Mobile Procedure Light
255-006	Wall Mount

### Accessories

002-1865-00 <sup>2,3</sup>	10' ceiling mount adapter kit
9A584001 <sup>1</sup>	Conversion kit w/cover (single to double light)

<sup>1</sup> - To convert an existing single light to a dual light, this conversion kit must be ordered in addition to a second 8' or 9' light

<sup>2</sup> - Converts an 8' or 9' light to a 10' ceiling height

<sup>3</sup> - One required per light arm

The following model numbers are recommended based on the following ceiling heights (measured from the bottom of the ceiling surface to the floor):

Less than 8'	255-005	Mobile Procedure Light
8' to 9'	255-001, 255-003	8' Single or Dual Procedure Light
9' to 10'	255-002, 255-004	9' Single or Dual Procedure Light
Greater than 10'	255-001, 255-003, 255-005	8' Single or Dual Light with 10' ceiling mount adapter kit(s) or Mobile Procedure Light



# Ritter 253

## LED Exam Light

- LED technology for up to 50,000 hours of useful life
- Designed to consume 45-60% less energy than other lights
- 4,500 K color temperature
- 1,700 fc intensity (18,298 lux) at 16-inch focal length
- Sleek, ergonomic and easy-to-clean design
- Mist color weight of light unit: 6 lb (253-002 light kit including bracket)
- Shipping weight: 8 lb (253-002 light kit including bracket)
- UL6060-1 and CAN/CSA C22.2 No. 601.1-M90



Exam Table Kit

### Mobile Caster Stand Kit

Model	Description
253-011	LED Exam Light - light only - excludes mounting hardware
253-012 <sup>3</sup>	LED Exam Light - exam table kit - models 204/22X*/62X* - includes light, splines and table mounting hardware, field installed
253-013 <sup>2</sup>	LED Exam Light - exam table kit - model 626 for 32" upholstery and Patient Support Rails - includes light, splines and table mounting hardware, field installed
253-015	LED Exam Light - procedure chair kit - models 230/641/646/647 - includes light, base rail and chair mounting hardware, field installed
253-016	LED Exam Light - mobile caster stand kit - includes light and caster stand
253-017	LED Exam Light - wall mount kit - includes light, wall mounting bracket and hardware, field installed

### Miscellaneous

9A619001	LED exam table mounting hardware (only) - models 204/22X*/604/62X* - field installed
9A620001	LED exam table mounting hardware (only) - model 626 w/ 32" upholstery - field installed
9A622001	LED procedure chair mounting hardware (only) - models 230/641/646/647 - field installed
9A624001	LED 630 procedure table spline mounting hardware (only), 28" wide upholstery, field installed
9A625001	LED 630 procedure table spline mounting hardware (only), 32" wide upholstery, field installed
9A626001 <sup>1</sup>	LED 630 procedure table rail mount hardware (only), field installed
9A459001	Disposable Sterility Sleeve; 5 boxes of 25 per carton
9A628001	LED light mobile caster stand (only)
9A623001	LED wall mount bracket (only)

\* Mounts to tables with a 2006 or newer manufacture date

<sup>1</sup> - Requires (1) accessory rails 9A557001 or 9A557002

<sup>2</sup> - Use 253-013 with Patient Support Rails™

<sup>3</sup> - Slight interference with Patient Support Rails, use 253-013 to minimize



# Ritter 250

## LED Exam Light

- LED technology for up to 25,000 hours of useful life
- Designed to consume 45-60% less energy than other lights
- 4,000 K color temperature
- 1,800 fc intensity (19,375 lux)
- UL6060-1 and CAN/CSA C22.2 No. 601.1-M90
- Mist color exam light weight: 20 lb
- Shipping weight: 25 lb



Caster Base Accessory

Model	Description
250-003	LED exam light
9A621001	Caster base accessory

# Midmark Diagnostic Medical Products (ECG, Holter, Spirometer, Stress, and Vital Signs) Limited Warranty

## SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will repair or replace devices manufactured by Midmark (except for those not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of applicable devices. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which, upon examination by Midmark, prove to be defective. This warranty extends only to the original retail purchaser of a product, and is not transferable or assignable.

Midmark warrants to the original retail purchaser that it will repair or replace software contained within the products manufactured by Midmark (except for those not warranted under "Exclusions") if: (1) the media on which the software is furnished exhibits defects in material or workmanship under normal use; or (2) the software does not substantially conform to its published specifications. Midmark does not warrant that the software: (1) is error free; (2) can be used without problems or interruptions; or (3) is free from vulnerability to intrusion or attack by viruses or other methods. Midmark also agrees to provide the original retail purchaser with software updates, if any, that become generally available to the public within the warranty period stated, for the specific product owned by the purchaser.

## APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original retail purchaser, shall be as follows:

- (a) One (1) year for software and diagnostic devices (unless otherwise noted below and except items listed under "Exclusions"),
- (b) Two (2) years for vital signs devices,
- (c) One (1) year for reusable blood pressure cuffs and SpO<sub>2</sub> sensors,
- (d) Ninety (90) days for all other accessories (sold in a product kit or purchased separately),

## OBTAINING WARRANTY SERVICE

Warranty service must be obtained through Midmark. Contact Midmark Technical Service for warranty service inquiries or issues by phone at 1-800-624-8950, option 2, or via email at techsupport@midmark.com.

It is the original retail purchaser's obligation to arrange for delivery of a product to Midmark for warranty service or exchange. It is also the original retail purchaser's obligation to comply with the warranty service instructions provided by Midmark.

## EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for the following:

- (1) any product that has been repaired, modified, or tampered with by any person or entity other than Midmark;
- (2) any product for which any label or tag attached to the product has been removed or tampered with;
- (3) defects, damage, or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, tampering, or failure to seek and obtain repair or replacement in a timely manner;
- (4) products which are not installed, used, and properly cleaned and maintained as stated in the Midmark Installation, User Guide and/or Operation Manual for the applicable product;
- (5) products considered to be of a consumable nature;
- (6) accessories or parts covered by the respective original manufacturers' warranty, including but not limited to treadmills, computers, printers, scales, carts, mounting options, and carrying cases, after the initial 90 days;
- (7) charges by anyone for adjustments, repairs, replacement parts, installation, or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (8) costs and expenses of routine maintenance and cleaning;
- (9) representations and warranties made by any person or entity other than Midmark; and
- (10) items provided through sales promotions.

## EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF DATA, LOSS OF USE, DOWNTIME, COVER, AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS, AND BENEFITS.

## NO AUTHORIZATION

No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

## WARRANTY DISCLAIMER

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

THIS WARRANTY WILL BE VOID IF ANY COVERED PRODUCT IS REPAIRED, MODIFIED, OR TAMPERED WITH BY ANY PERSON OR ENTITY OTHER THAN MIDMARK, OR IF ANY LABEL OR TAG ATTACHED TO A COVERED PRODUCT IS REMOVED OR TEMPERED WITH.

## STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period. In the event multiple warranty periods exist with respect to a product, the ninety (90) day period provided for herein shall begin to run from expiration of the warranty period for the component to which the claim relates.

99-96-00750014 Rev G

## Midmark Digital Diagnostic Devices Part Number Ordering Information

Note: Midmark recommends that all Midmark devices drop ship directly to the customer

### TO ORDER:

Midmark devices require either an electronic medical record (EMR) that has integrated Midmark software or Midmark IQmanager® Software to operate. Devices include, but are not limited to, the part numbers shown below:

4-000-0025  
4-000-0026  
4-000-0062  
4-000-0110  
4-000-0113  
4-000-0116

If the customer *currently has* an EMR software package that has integrated Midmark Software, no additional software is required.

Place an order for the appropriate device(s) and please provide:

- EMR name and version
- Customer facility + contact name, address and phone
- Customer computer environment (Thin Client, Fat Client)

If the customer *does not have* an EMR, or **is** using Allscripts TouchWorks™ EHR, Epic® EMR or GE Centricity®, please contact your local Midmark Sales Representative to order the appropriate Midmark software solution. Additional charges or license fees will apply.

For ordering questions regarding Midmark devices, please contact your local Midmark Sales Representative or call 800.624.8950, option 2.

Allscripts TouchWorks EHR, Epic EMR and GE Centricity are not trademarks of Midmark Corporation and are the sole property of their respective owners.



## Midmark HEINE Diagnostic Instruments Limited Warranty

### SCOPE OF WARRANTY

Midmark Corporation ("Midmark") warrants to the original retail purchaser that it will, at Midmark's option, repair or replace diagnostic instruments manufactured by HEINE and sold by Midmark (except for those not warranted under "Exclusions") that are defective in material or workmanship under normal use and service. The sole remedy under this limited warranty is the repair or replacement, at Midmark's option, of the applicable components. This limited warranty shall only apply to defects that are reported to Midmark within the applicable warranty period and which are determined to exist upon examination by Midmark. This warranty extends only to the original retail purchaser of a product and is not transferable or assignable. Replacement components or products may be used and/or refurbished components or products, provided they are of like quality and specifications as new components or products.

### APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of invoice to the original retail purchaser of the product shall be as follows:

- (1) Ten (10) years for Diagnostic instruments (other than GAMMA G5 Sphygmomanometers) manufactured by HEINE and sold by Midmark;
- (2) Five (5) years for GAMMA G5 Sphygmomanometers manufactured by HEINE and sold by Midmark; and
- (3) Twelve (12) months on parts and labor.

### OBTAINING WARRANTY SERVICE

Warranty service must be obtained through Midmark. Midmark may be contacted for warranty service inquiries or issues via email at [www.midmark.com](http://www.midmark.com); by phone at 1-800-MIDMARK or by mail to Midmark Corporation, 60 Vista Drive, Versailles, Ohio 45380. It is the retail purchaser's obligation to arrange for delivery of a product to Midmark or one of its authorized dealers for warranty service, which delivery shall be at retail purchaser's expense. It is also the retail purchaser's obligation to comply with the warranty service instructions provided either by Midmark or its authorized dealer. The retail purchaser must provide Midmark with completed warranty registration information within thirty (30) days after purchase in order to obtain the benefits of this limited warranty.

### EXCLUSIONS:

This limited warranty does not cover and Midmark shall not be liable for the following:

- (1) Defects, damage or other conditions caused, in whole or in part, by misuse, abuse, negligence, alteration, accident, freight damage, negligent storage, tampering or failure to seek and obtain repair or replacement in a timely manner;
- (2) Products which are not installed, used, and properly cleaned and maintained as required or recommended in the Midmark "Installation" and/or "Installation/Operation Manual" for the applicable product, including the specified structural and operational environmental conditions and electrical requirements;
- (3) Normal wear and tear
- (4) Fire or theft
- (5) Products considered to be single-use items or consumables (e.g. bulbs, cuffs, and batteries)
- (6) Rechargeable batteries and batteries
- (7) Use of non-original HEINE spare parts
- (8) Exposed and accessible control and supply cables with corresponding plug connection
- (9) Charges by anyone for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which are not expressly authorized in writing in advance by Midmark;
- (10) Costs and expenses of routine maintenance and cleaning;
- (11) Representations and warranties made by any person or entity other than Midmark.
- (12) Custom manufactured products;
- (13) Alterations or modifications to the product by any person or entity other than Midmark; and
- (14) Products that would otherwise be covered under this limited warranty, but are acquired: (i) from a person or entity that is not Midmark or one of its authorized dealers; or (ii) from a Midmark dealer that is not authorized to sell the product at issue in the geographic territory where the purchaser is located, or is not authorized to sell the product at issue within the medical market

### EXCLUSIVE REMEDY; CONSEQUENTIAL DAMAGES DISCLAIMER

MIDMARK'S ONLY OBLIGATION UNDER THIS LIMITED WARRANTY IS THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS. MIDMARK SHALL NOT BE LIABLE FOR AND HEREBY DISCLAIMS ANY DIRECT, SPECIAL, INDIRECT, INCIDENTAL, EXEMPLARY OR CONSEQUENTIAL DAMAGES OR DELAYS, INCLUDING, BUT NOT LIMITED TO, DAMAGES FOR LOSS OF PROFITS OR INCOME, LOSS OF USE, LOSS OF DATA, DOWNTIME, COVER AND EMPLOYEE OR INDEPENDENT CONTRACTOR WAGES, PAYMENTS AND BENEFITS. THIS DISCLAIMER SHALL SURVIVE ANY FAILURE OR ASSERTED FAILURE OF THE ESSENTIAL PURPOSE OF THIS LIMITED WARRANTY OR ITS REMEDIES SPECIFIED HEREIN.

### WARRANTY DISCLAIMER

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

### STATUTE OF LIMITATIONS

No action may be brought against Midmark for breach of this limited warranty, an implied warranty, if any, or for any other claim arising out of or relating to the products, more than ninety (90) days following expiration of the limited warranty period.

### NO AUTHORIZATION

No person or firm is authorized to create or approve for Midmark any other obligation or liability in connection with the products.

\*Warranty information included in product manuals

## Midmark Medical Diagnostic Products Support Information

For technical assistance with a Midmark medical diagnostic product, please contact Technical Service at 844.856.1230.

Select: Option 2 – for IQvitals® and IQvitals® Zone™ Digital Vital Signs, IQspiro® Digital Spirometer, IQecg® Digital ECG, IQholter® Digital Holter, Software, Supplies and Midmark IQcart® Diagnostic Equipment Cart Available Monday-Friday, 8:00 am – 8:00 pm, EST.

Select: Option 3 – for Diagnostic Instrument Wall/Desktop Kits, Instrument Replacement Parts and Supplies (Otoscope/Ophthalmoscope) and Sphygmomanometer Kits. Available 8:00 am – 7:00 pm Monday-Thursday and 8:00 am – 5:00 pm Friday, EST.

For troubleshooting tips and product use information, refer to the operation manual. Additional help topics, self-guided user tools and documentation can be accessed at any time through our support section and knowledge base at [midmark.com](http://midmark.com).

### Return Material Authorization

To return any product for repair, an RMA (Return Material Authorization) number must be obtained from a Midmark Technical Service Specialist. The RMA number should be referenced on the package(s) containing the items to be returned.

### Shipping

Before shipping any product to Midmark, be certain that an RMA number has been issued and that all guidelines regarding this authorization are followed. We highly recommend that you follow all guidelines for the shipment of medical products set forth by the shipping company you choose to use. If a question should arise regarding the appropriate method of shipment, please feel free to ask when calling for your RMA number. It is the responsibility of the customer when shipping a product to ensure that all packages and their contents get to Midmark safely.

**Midmark will not assume responsibility for damage incurred by improper packaging, shipment or product use. Such actions will void all applicable warranties.**

# IQvitals® Zone™ Vital Signs Wireless Monitor

with Touchscreen, Serenity BP™ and Exergen TemporalScanner™

- includes Spot, Averaging, Continuous, Interval and SPRINT Protocol BP Modes



Shown with IQvitals Zone Vital Signs Monitor mobile cart option\*

\*Mobile cart not included in pricing.

- 1 - IQvitals Zone Wireless Monitor
- 1 - Reusable Dual Tube BP Cuff – Adult (26-35cm)
- 1 - Reusable Dual Tube BP Cuff – Large Adult (32-42cm)
- 1 - Dual Tube BP Hose – 6.5'
- 1 - Exergen TemporalScanner Thermometer
- 1 - Exergen Disposable Caps (10 ea)
- 1 - IQvitals Zone Wireless USB Dongle
- 1 - IQvitals Zone USB Cable – 10'
- 1 - IQvitals Zone Laptop Tray – 15"
- 1 - Front Storage Bin Kit
- 1 - Device to Temperature Module Cable
- 1 - Three-Arm Knobs (Set of 2)
- 1 - AC Power Supply
- 1 - Setup Guide
- 1 - Quick Reference Guide
- 1 - Operation Manual CD

Model	Description
4-000-0724 <sup>1</sup>	IQvitals Zone Wireless Monitor, with Serenity BP, Exergen TemporalScanner and Masimo® SpO <sub>2</sub>
4-000-0722 <sup>1</sup>	IQvitals Zone Wireless Monitor, with Serenity BP, Exergen TemporalScanner and Nellcor™ SpO <sub>2</sub>
4-000-0720 <sup>1</sup>	IQvitals Zone Wireless Monitor, with Serenity BP and Exergen TemporalScanner

<sup>1</sup>Mounting option is required. See available options below.

## IQvitals Zone Monitor Mounting Options

Model	Description
2-200-0072	IQvitals Zone Mobile Cart
2-200-0082	IQvitals Zone Wall Mount Articulating Arm

## Midmark Medical Diagnostic Products Support Information

For technical assistance with a Midmark medical diagnostic product, please contact Technical Service at 844.856.1230.

Select: Option 2 – for IQvitals® and IQvitals® Zone™ Digital Vital Signs, IQspiro® Digital Spirometer, IQecg® Digital ECG, IQholter® Digital Holter, Software, Supplies and Midmark IQcart® Diagnostic Equipment Cart Available Monday-Friday, 8:00 am – 8:00 pm, EST.

Select: Option 3 – for Diagnostic Instrument Wall/Desktop Kits, Instrument Replacement Parts and Supplies (Otoscope/Ophthalmoscope) and Sphygmomanometer Kits. Available 8:00 am – 7:00 pm Monday-Thursday and 8:00 am – 5:00 pm Friday, EST.

For troubleshooting tips and product use information, refer to the operation manual. Additional help topics, self-guided user tools and documentation can be accessed at any time through our support section and knowledge base at [midmark.com](http://midmark.com).

### Return Material Authorization

To return any product for repair, an RMA (Return Material Authorization) number must be obtained from a Midmark Technical Service Specialist. The RMA number should be referenced on the package(s) containing the items to be returned.

### Shipping

Before shipping any product to Midmark, be certain that an RMA number has been issued and that all guidelines regarding this authorization are followed. We highly recommend that you follow all guidelines for the shipment of medical products set forth by the shipping company you choose to use. If a question should arise regarding the appropriate method of shipment, please feel free to ask when calling for your RMA number. It is the responsibility of the customer when shipping a product to ensure that all packages and their contents get to Midmark safely.

**Midmark will not assume responsibility for damage incurred by improper packaging, shipment or product use. Such actions will void all applicable warranties.**

# IQvitals® Zone™ Vital Signs Wireless Monitor

with Touchscreen, Serenity BP™ and Exergen TemporalScanner™

- includes Spot, Averaging, Continuous, Interval and SPRINT Protocol BP Modes



Shown with IQvitals Zone Vital Signs Monitor mobile cart option\*

\*Mobile cart not included in pricing.

- 1 - IQvitals Zone Wireless Monitor
- 1 - Reusable Dual Tube BP Cuff – Adult (26-35cm)
- 1 - Reusable Dual Tube BP Cuff – Large Adult (32-42cm)
- 1 - Dual Tube BP Hose – 6.5'
- 1 - Exergen TemporalScanner Thermometer
- 1 - Exergen Disposable Caps (10 ea)
- 1 - IQvitals Zone Wireless USB Dongle
- 1 - IQvitals Zone USB Cable – 10'
- 1 - IQvitals Zone Laptop Tray – 15"
- 1 - Front Storage Bin Kit
- 1 - Device to Temperature Module Cable
- 1 - Three-Arm Knobs (Set of 2)
- 1 - AC Power Supply
- 1 - Setup Guide
- 1 - Quick Reference Guide
- 1 - Operation Manual CD

Model	Description
4-000-0724 <sup>1</sup>	IQvitals Zone Wireless Monitor, with Serenity BP, Exergen TemporalScanner and Masimo® SpO <sub>2</sub>
4-000-0722 <sup>1</sup>	IQvitals Zone Wireless Monitor, with Serenity BP, Exergen TemporalScanner and Nellcor™ SpO <sub>2</sub>
4-000-0720 <sup>1</sup>	IQvitals Zone Wireless Monitor, with Serenity BP and Exergen TemporalScanner

<sup>1</sup>Mounting option is required. See available options below.

## IQvitals Zone Monitor Mounting Options

Model	Description
2-200-0072	IQvitals Zone Mobile Cart
2-200-0082	IQvitals Zone Wall Mount Articulating Arm

# IQvitals® Vital Signs Device

with Touchscreen - Blood Pressure and Temperature



- 1 - IQvitals device w/Blood Pressure and Temperature
- 1 - AC Power Supply
- 1 - Reusable Blood Pressure Cuff - Adult
- 1 - Reusable Blood Pressure Cuff - Large Adult
- 1 - Blood Pressure Hose - 6.5'
- 1 - Oral/Axillary Temperature Probe
- 1 - Temperature Probe Covers (box of 20)
- 1 - USB Cable - 10'
- 1 - Operation Manual CD
- 1 - Quick Reference Guide

Model	Description
4-000-0500*	IQvitals device with Blood Pressure and Alaris® Turbo Temp® Thermometer

## IQvitals Vital Signs Device with Touchscreen - Blood Pressure, Temperature and SpO<sub>2</sub>

- 1 - IQvitals device w/Blood Pressure, Temperature and SpO<sub>2</sub>
- 1 - AC Power Supply
- 1 - Reusable Blood Pressure Cuff - Adult
- 1 - Reusable Blood Pressure Cuff - Large Adult
- 1 - Blood Pressure Hose - 6.5'
- 1 - Reusable SpO<sub>2</sub> Adult Finger Sensor
- 1 - SpO<sub>2</sub> Extender Cable - 4'
- 1 - Oral/Axillary Temperature Probe
- 1 - Temperature Probe Covers (box of 20)
- 1 - USB Cable - 10'
- 1 - Operation Manual CD
- 1 - Quick Reference Guide

Model	Description
4-000-0510*	IQvitals device with Blood Pressure, Alaris Turbo Temp Thermometer and SpO <sub>2</sub>

\*Available to the U.S. market only

# IQvitals® Accessories and Supplies

Model	Description
1-100-1603	Fairbanks® TeleWeigh® Digital Scale
1-100-1605	IQvitals Vital Signs Thermal Printer
3-004-2004	IQvitals Vital Signs Mobile Cart Printer Bracket
3-004-2010	IQvitals Vital Signs Mobile Cart Scale Mount
3-009-0010	AC Power Supply
3-009-0012	Power Supply Adapter Kit
3-009-0014	IQvitals + IQvitals PC Replacement Lithium ion battery
3-009-0016	USB Cable - 10'
3-009-0018	Serial Cable - 10'
3-009-0022	Blood Pressure Hose - 6.5'
3-009-0100	Blood Pressure Hose - 10'
3-009-0021	Reusable SpO <sub>2</sub> Small Finger Sensor, Envitec®
3-009-0020	Reusable SpO <sub>2</sub> Adult Finger Sensor, Envitec
3-009-0026	SpO <sub>2</sub> Extender Cable, 4', Envitec
3-009-0024*	Alaris® Turbo Temp® Oral/Axillary Temperature Probe
3-009-0058*	Turbo Temp Probe Covers, Carton of 10 boxes (20 covers per box)
3-009-0068	Reusable Blood Pressure Cuff, Infant (8-14 cm)
3-009-0070	Reusable Blood Pressure Cuff, Child (13-20 cm)
3-009-0062	Reusable Blood Pressure Cuff, Small Adult (18-26 cm)
3-009-0064	Reusable Blood Pressure Cuff, Adult, Midmark (26-35 cm)
3-009-0066	Reusable Blood Pressure Cuff, Large Adult, Midmark (32-42 cm)
3-009-0072	Reusable Blood Pressure Cuff, Adult Long (29-38 cm)
3-009-0074	Reusable Blood Pressure Cuff, Large Adult Long (35-44 cm)
3-009-0076	Reusable Blood Pressure Cuff, Thigh (42-50 cm)
3-004-2000	IQvitals Vital Signs Mobile Cart
3-009-0001	IQvitals Vital Signs Countertop Mount
3-009-0003	IQvitals Vital Signs Wall Mount with Basket
3-004-2008	IQvitals Vital Signs Equipment Pole Mount

\*Available to the U.S. market only

## Connectivity Supplies for Midmark 626

Model	Description
9A551002	IQvitals Vital Signs Table Mount, non-hub
9A478001	IQscale serial cable for IQvitals Vital Signs, 6' length, coiled
9A478002	IQscale serial cable for IQvitals Vital Signs, 15' length, straight
9A478003	IQscale serial cable for IQvitals Vital Signs, 30' length, straight
9A478004	IQscale serial cable for IQvitals Vital Signs, 50' length, straight
9A629001	IQvitals Vital Signs table mount with USB and serial cables
9A629002	IQvitals Vital Signs table mount with serial cable
9A629003	IQvitals Vital Signs table mount, non-hub
002-10064-00	626 Chair Connectivity Kit*

\*Not recommended for connecting Midmark Digital Vital Signs to PC running Windows® 10 64-bit OS



# Midmark IQspiro® Digital Spirometer

- 1 - IQspiro device
- 10 - Disposable Spirometer Mouthpieces
- 1 - 10/pack Disposable Nose Clips
- 1 - Mouse Pad
- 1 - Operation Manual CD
- 1 - Quick Reference Guide
- 1 - Carrying Case

Model	Description
4-000-0026	IQspiro - USB version with syringe
4-000-0025	IQspiro - USB version

---

## Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

---

## IQspiro Accessories and Supplies

Model	Description
1-100-0007	3-Liter Spirometer Calibration Syringe
2-100-1205	Disposable Spirometer Mouthpieces (box of 25)
2-100-1206	Disposable Spirometer Mouthpieces (case - 4 boxes)
2-100-1207	Box of 500 Mouthpieces, individually bagged, bulk
2-244-0001	Disposable Nose Clips (box of 25)





# Midmark IQecg® Digital ECG

## with Lead Management (Interpretive)

- 1 - IQecg Acquisition Module
- 1 - 10-Lead Patient Cable with Lead Management
- 1 - 10/pack Clear ECG Clips
- 1 - 100/pack Disposable EZ-Trode ECG Electrodes
- 1 - Mouse Pad
- 1 - Operation Manual CD
- 1 - Quick Reference Guide
- 1 - Carrying Case

Model	Description
4-000-0062	IQecg device - with Lead Management

## Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

## IQecg Supplies

Model	Description
3-047-0001*	Universal ECG Clips (3mm, 4mm + Snap) 10/pack
3-047-0005	Clear ECG Clips (4mm + Snap) 10/pack
3-100-0203	Patient Cable, with Lead Management, IQecg
2-100-0205	Disposable EZ-Trode ECG Electrodes (box of 500)
2-100-0206	Disposable EZ-Trode ECG Electrodes (case - 10 boxes)

\*Available to the U.S. market only

# Midmark IQholter® Digital Holter



## Holter Kit with Digital Recorder

- 1 - IQholter Recorder
  - Includes: lead wire set, secure digital card
- 1 - Reusable IQholter Recorder Pouch and Belt
- 1 - Secure Digital Card Reader
- 1 - Software Security Key, USB Version
- 4 - Holter Prep Kits
- 1 - Mouse Pad
- 1 - Operation Manual CD
- 1 - Quick Reference Guide

Model	Description
4-000-0110*	IQholter with one (1) Recorder
4-000-0113*	IQholter EX with one (1) Recorder
4-000-0116*	IQholter EP with Pacemaker detection and one (1) Recorder

\*Available to the U.S. market only

## Midmark Digital Holter Recorder

- 1 - IQholter Recorder
- 1 - Lead wire set
- 1 - Secure Digital Card
- 1 - Operation Manual CD
- 4 - Holter Prep Kits
- 1 - Reusable Digital Holter Recorder Pouch and Belt

Model	Description
4-000-0119*	IQholter Recorder
1-370-0020	IQholter Recorder Upgrade Kit (24, 48, 72 Hours) One Upgrade Kit Required per Recorder

## Software

Please contact your local Midmark Sales Representative to order the appropriate Midmark software solution, including the new IQmanager® version 10.0.0 for non-EMR users.

## Secure Memory Card

Model	Description
1-370-0010	IQholter® Secure Digital Memory Card
1-370-0004	USB Secure Digital Card Reader

## IQholter Supplies

Model	Description
2-100-0080	Prep Kit - Digital Holter, 5-lead (each) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 abrading pad, 2 alcohol wipes, 5 electrodes
2-100-0081	Prep Kit - Digital Holter, 5-lead (case) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 abrading pad, 2 alcohol wipes, 5 electrodes (case - 30 kits)
2-100-0090	Prep Kit - Digital Holter, 7-lead (each) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 abrading pad, 2 alcohol wipes, 7 electrodes
2-100-0091	Prep Kit - Digital Holter, 7-lead (case) Each kit contains: 1 diary, 2 AA batteries, 1 razor 1 abrading pad, 2 alcohol wipes, 7 electrodes (case - 30 kits)
90-30-0205*	Lead Wire, 5-lead, Midmark IQholter (Blue)
90-30-0206*	Lead Wire, 7-lead, Midmark IQholter (Black)
2-100-0060	Midmark Holter Electrodes (box of 50)
2-100-0050	Disposable IQholter Recorder Pouches (pack of 25)
2-100-0051	Reusable IQholter Recorder Pouch with Belt

\*Available to the U.S. market only



## Midmark HEINE Diagnostic Instrument Wall Kits

- 1 - EN200® Wall Transformer
- 1 - Otoscope with polycarbonate, 11.5' (3.5 m) corded handle
- 1 - Ophthalmoscope with polycarbonate, 11.5' (3.5 m) corded handle

Model	Description
A-095-12-208-166	EN200 Wall Transformer with BETA® 400 LED Otoscope and BETA® 200 LED Ophthalmoscope
A-095-12-209-166	EN200 Wall Transformer with BETA 200 LED Otoscope and BETA 200 LED Ophthalmoscope
A-095-12-204-166	EN200 Wall Transformer with BETA® 400 XHL Otoscope and BETA® 200 XHL Ophthalmoscope
A-095-12-205-166	EN200 Wall Transformer with BETA 200 XHL Otoscope and BETA 200 XHL Ophthalmoscope
A-095-12-210-166	EN200 Wall Transformer with K180® XHL Otoscope and K180® XHL Ophthalmoscope

---

## Midmark HEINE Diagnostic Instrument Wall Kit Accessories

Model	Description
X-095-12-006-166	Wallboard (for device combination of EN200 + Tip Dispenser - no instruments)
X-095-12-007-166	Wallboard (for device combination of EN200 + Tip Dispenser + GAMMA® XXL Sphygmomanometer - no instruments)
B-000-11-149-166	AllSpec® Tip Dispenser - white



## Midmark HEINE Diagnostic Instrument Desktop Kits

- 1 – NT4 Desktop Charger
- 1 – Otoscope with BETA®4 NT, Lithium-ion rechargeable handle
- 1 – Ophthalmoscope with BETA®4 NT, Lithium-ion rechargeable handle

Model	Description
A-853-24-420-166	NT4 Desktop Charger with BETA® 400 LED Otoscope and BETA® 200 LED Ophthalmoscope
A-832-24-420-166	NT4 Desktop Charger with BETA 200 LED Otoscope and BETA 200 LED Ophthalmoscope
A-853-23-420-166	NT4 Desktop Charger with BETA® 400 XHL Otoscope and BETA® 200 XHL Ophthalmoscope
A-832-23-420-166	NT4 Desktop Charger with BETA 200 XHL Otoscope and BETA 200 XHL Ophthalmoscope
A-879-23-420-166	NT4 Desktop Charger with K180® XHL Otoscope and K180 XHL Ophthalmoscope

## Midmark HEINE Diagnostic Instrument Desktop Kit Accessories

Model	Description
B-000-11-149-166	AllSpec® Tip Dispenser - white



## Midmark HEINE Diagnostic Instrument Replacement Parts

Model	Description
B-008-11-400-166	BETA® 400 LED Otoscope (instrument head only)
B-002-11-400-166	BETA® 400 XHL Otoscope (instrument head only)
B-008-11-500-166	BETA® 200 LED Otoscope (instrument head only)
B-002-11-500-166	BETA® 200 XHL Otoscope (instrument head only)
B-002-11-550-166	K180® XHL Otoscope (instrument head only)
C-008-30-100-166	BETA 200 LED Ophthalmoscope (instrument head only)
C-002-30-100-166	BETA 200 XHL® Ophthalmoscope (instrument head only)
C-002-30-207-166	K180 XHL Ophthalmoscope (instrument head only)
X-002-99-494-166	NT4 Desktop Charger (without instruments)
X-007-99-396-166	BETA® 4 NT Rechargeable Handle (with battery and handle bottom cap)
X-002-99-394-166	BETA® 4 NT Handle Bottom Cap
X-095-12-250-166	EN200® Wall Transformer (without instruments)



## Midmark HEINE Diagnostic Instrument Supplies

Model	Description
B-000-11-127-166	AllSpec® Disposable Tips 4 mm (box of 1,000)
B-000-11-137-166	AllSpec® Disposable Tips 4 mm (case of 10 boxes, 1,000/box)
B-000-11-128-166	AllSpec® Disposable Tips 2.5 mm (box of 1,000)
B-000-11-138-166	AllSpec® Disposable Tips 2.5 mm (case of 10 boxes, 1,000/box)
X-002-88-078-166	XHL® Xenon Halogen 3.5 V Bulb for BETA 400, BETA 200 and K180 Otoscopes (package of 6)
X-002-88-070-166	XHL Xenon Halogen 3.5 V Bulb for BETA 200 Ophthalmoscope (package of 6)
X-002-88-086-166	XHL Xenon Halogen 3.5 V Bulb for K180 Ophthalmoscope (package of 6)
X-007-99-383-166	Rechargeable Lithium-ion Battery



## Midmark HEINE GAMMA<sup>®</sup> XXL LF Sphygmomanometer Kit

- 1 – GAMMA XXL LF Professional Sphygmomanometer
- 1 – Twin-tube adult cuff (29 – 41 cm circ.)
- 1 – Twin-tube adult large cuff (35 – 47 cm circ.)
- 1 - Universal mounting plate

Model	Description
M-000-09-323-166	GAMMA XXL LF Sphygmomanometer with universal mounting plate and two twin-tube cuffs (adult and adult large)

## Midmark HEINE GAMMA XXL Accessories

Model	Description
M-000-09-701-166	Twin-tube LF blood pressure cuff, Infant (6 x 28 cm; 8 – 13 cm circumference)
M-000-09-703-166	Twin-tube LF blood pressure cuff, Child (10 x 34 cm; 13 - 20 cm circumference)
M-000-09-705-166	Twin-tube LF blood pressure cuff, Adult Small (14 x 53 cm; 20 – 29 cm circumference)
M-000-09-629-166	Twin-tube LF blood pressure cuff, Adult (14 x 58 cm; 29 – 41 cm circumference)
M-000-09-631-166	Twin-tube LF blood pressure cuff, Adult Large (17 x 70 cm; 35 - 47 cm circumference)
M-000-09-625-166	Twin-tube LF blood pressure cuff, Thigh (20 x 86 cm; >41 cm circumference)



## Midmark HEINE GAMMA® G5 Sphygmomanometer Kit

- 1 – GAMMA G5 Standard Sphygmomanometer
- 1 – Single-tube adult cuff (29 – 41 cm circ.)
- 1 – Single-tube adult large cuff (35 – 47 cm circ.)

Model	Description
M-000-09-230-166	GAMMA G5 Sphygmomanometer with two single-tube cuffs (adult and adult large)

---

## Midmark HEINE GAMMA G5 Accessories

Model	Description
M-000-09-700-166	Single-tube LF blood pressure cuff, Infant (6 x 28 cm; 8 – 13 cm circumference)
M-000-09-702-166	Single-tube LF blood pressure cuff, Child (10 x 34 cm; 13 - 20 cm circumference)
M-000-09-704-166	Single-tube LF blood pressure cuff, Adult Small (14 x 53 cm; 20 – 29 cm circumference)
M-000-09-619-166	Single-tube LF blood pressure cuff, Adult (14 x 58 cm; 29 – 41 cm circumference)
M-000-09-706-166	Single-tube LF blood pressure cuff, Adult Large (17 x 70 cm; 35 - 47 cm circumference)
M-000-09-615-166	Single-tube LF blood pressure cuff, Thigh (20 x 86 cm; >41 cm circumference)





# Midmark IQcart® Diagnostic Equipment Cart

(Some assembly required)

- 1 - Storage Basket
- 1 - IQecg® Arm
- 1 - Digital Spirometer Holster
- 1 - Set of Locking Disks

Model	Description
3-004-1000	Midmark IQcart with Arm

## IQcart Accessories and Supplies

Model	Description
3-004-2008	IQvitals® Equipment Pole Mount
3-004-6028	IQcart Basket (8" x 6" x 8")

## Midmark Workstations: Limited Warranty

### SCOPE OF WARRANTY

Midmark Corporation warrants to the original purchaser its Midmark Workstation products manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the parts or the products the defects of which are reported to Midmark within the applicable warranty period and which upon examination by Midmark prove to be defective.

### APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

- A. Five years for all Midmark Workstation products and components (except for the shorter periods in B), including drawer fronts and casters
- B. One year for Power system components and coiled power cords

### EXCLUSIONS

This warranty does not cover and Midmark shall not be liable for: (1) Failure to install or use according to the Midmark user and installation manuals (2) Normal wear and tear, repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage or tampering (including, without limitation, use of product in unsuitable environments or conditions); (3) matching of color, grain or texture except to commercially acceptable standards; (4) changes in color caused by natural or artificial light; (5) products which are

not installed, used or properly cleaned as required in Midmark's written Installation Manual; (6) specially manufactured products; (7) products considered to be of a consumable nature such as light bulbs and surge suppression products; (8) accessories or parts not manufactured by Midmark or use of non-authorized Midmark components; (9) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark. 10) Care Exchange products and accessories (pgs 69-83) follow warranty on page 5.

### EXCLUSIVE REMEDY

Midmark's only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay.

### NO AUTHORIZATION

No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.



## Midmark 6211 Secure Laptop Workstation

- 26" of height adjustment from 23" to 49"
- Keyed laptop security system
- 23" D x 26" W worksurface
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- 19.5" D x 20.5" W base dimension
- Weight capacity: 20 lb
- Unit weight: 58 lb
- Shipping weight: 103 lb

**Model**  
6211-001-xxx

**Description**  
Secure Laptop Workstation

<sup>1</sup> - See color code chart on page 71

### Accessories

Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/Steel Cable Mgmt-LH
PMSR	9A655001	Power Mount w/Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



## Midmark 6212 Flat Panel PC Workstation

- 26" of height adjustment from 23" to 49"
- 23" D x 26" W worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive tilt
- Monitor arm capacity range 5 lb min - 20 lb max
- Monitor arm independent rotation of 180°
- 22.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 55 lb
- Shipping weight: 100 lb

### Model

6212-001-xxx

### Description

Flat Panel PC Workstation

<sup>1</sup> - See color code chart on page 71

### Accessories

Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/Steel Cable Mgmt-LH
PMSR	9A655001	Power Mount w/Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Spiro Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



## Midmark 6213 Flat Panel Secure PC Workstation

- 26" of height adjustment from 23" to 49"
- Keyed PC security system
- 23" D x 26" W worksurface
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive tilt
- Monitor arm capacity range 5 lb min - 20 lb max
- Monitor arm independent rotation of 180°
- 22.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 64 lb
- Shipping weight: 109 lb

### Model

6213-001-xxx

### Description

Flat Panel Secure PC Workstation

<sup>1</sup> - See color code chart on page 71

### Accessories

Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PCML	9A669001	Plastic Bin - MG Coiled PWR Cord-LH
PCMR	9A669001	Plastic Bin - MG Coiled PWR Cord-RH
PMBL	9A655002	Power Mount w/Bin Cable Mgmnt-LH
PMBR	9A655002	Power Mount w/Bin Cable Mgmnt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHEP	3-004-6030	Spiro Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



## Midmark 6214 Procedure Workstation

- 26" of height adjustment from 23" to 49"
- 21" D x 26" W worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Includes: rear basket, BP cuff holder, integrated spiro holsters and equipment pole
- 19.5" D x 20.5" W base dimension
- Weight capacity: 30 lb
- Unit weight: 56 lb
- Shipping weight: 101 lb

**Model**  
6214-001-xxx

**Description**  
Procedure Workstation

<sup>1</sup> - See color code chart on page 71

### Accessories

Plant Installed	Field Installed	
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
IQEP	3-004-2008	Digital Vital Signs EP Mount
PMSL	9A655001	Power Mount w/Steel Cable Mgmt-LH
PMSR	9A655001	Power Mount w/Steel Cable Mgmt-RH
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHEP	3-004-6030	Spiro Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf



## Midmark 6215 Compact Workstation

- 26" of height adjustment from 23" to 49"
- 19.5" D x 20.5" W base dimension
- 21" D x 26" W worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 30 lb
- Unit weight: 48 lb
- Shipping weight: 93 lb

### Model

6215-001-xxx

### Description

Compact Workstation

<sup>1</sup> - See color code chart on page 71

### Accessories

Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLW	9A647001	Drawer w/Lock-Work Surf Mnt
DWNW	9A646001	Drawer-Work Surf Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
MPSF	9A648001	Mobile Printer Shelf
PMSL	9A655001	Power Mount w/Steel Cable Mgmt-LH
PMSR	9A655001	Power Mount w/Steel Cable Mgmt-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Spiro Holder EP Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



## Midmark 6217 Secure Laptop Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- 19.5" D x 20.5" W base dimension
- Keyed laptop security system
- 23" D x 22" W worksurface w/ slideout mouse surface
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Weight capacity: 30 lb
- Unit weight: 60 lb
- Shipping weight: 105 lb

**Model**  
6217-001-xxx

**Description**  
Secure Laptop Workstation + Mouse Tray

<sup>1</sup> - See color code chart on page 71

### Accessories

Plant Installed	Field Installed	
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DWLM	9A647001	Drawer w/Lock-Mouse Trk Mnt
DWNM	9A646001	Drawer-Mouse Track Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
MPSF	9A648001	Mobile Printer Shelf
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Spiro Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
USCB	9A656001	Adjustable Utility Shelf
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount





## Midmark 6218

### Flat Panel PC Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Monitor arm capacity range 5 lb min - 20 lb max
- Monitor arm independent rotation of 180°
- 23" D x 22" W worksurface w/ slideout mouse surface
- 22.5" D x 20.5" W base dimension
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. and up to 40° positive
- Weight capacity: 30 lb
- Unit weight: 55 lb
- Shipping weight: 100 lb

#### Model

6218-001-xxx

#### Description

Flat Panel PC Workstation + Mouse Tray

<sup>1</sup> - See color code chart on page 71

#### Accessories

Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLM	9A647001	Drawer w/Lock-Mouse Trk Mnt
DWNM	9A646001	Drawer-Mouse Track Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SCML	9A668001	Steel Brkt + MG Coiled Cord-LH
SCMR	9A668001	Steel Brkt + MG Coiled Cord-RH
SHEP	3-004-6030	Spiro Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



## Midmark 6219 Flat Panel Secure PC Workstation + Mouse Tray

- 26" of height adjustment from 23" to 49"
- Monitor arm capacity range 5 lb min - 20 lb max
- Keyed PC security system
- Monitor arm independent rotation of 180°
- 23" D x 22" W worksurface w/ slideout mouse surface
- 22.5" D x 20.5" W base dimension
- 11" D x 20.5" W secondary worksurface
- Two total-locking casters in the front lock both the wheel and swivel simultaneously (two non-locking casters in the rear)
- Monitor arm w/ 75 mm and 100 mm mounting patterns
- Monitor arm w/ 9.5" independent vertical adj. & up to 40° positive tilt
- Weight capacity: 30 lb
- Unit weight: 64 lb
- Shipping weight: 109 lb

### Model

6219-001-xxx

### Description

Flat Panel Secure PC Workstation + Mouse Tray

<sup>1</sup> - See color code chart on page 71

### Accessories

Plant Installed	Field Installed	
AIOM	N/A	Monitor over 10 lb
BCHD	9A664001	Bar Code Holder Mount
BCHL	9A659001	Bar Code Holder-LH
BCHR	9A659001	Bar Code Holder-RH
BPCH	9A658001	Blood Pressure Cuff Holder
CCWT	9A657001	Cart Counterweight
CPUB	9A663001	CPU Holder-Back
CPUF	9A663001	CPU Holder-Front
DBCM	9A653001	Double Basket
DMMT	9A667001	Dual Monitor Mount
DWLM	9A647001	Drawer w/Lock-Mouse Trk Mnt
DWNM	9A646001	Drawer-Mouse Track Mnt
EPWS	9A651001	Equipment Pole
IQEP	3-004-2008	Digital Vital Signs Equipment Pole Mount
LPDH	9A662001	Laptop Display Holder
MPSF	9A648001	Mobile Printer Shelf
PCML	9A669001	Plastic Bin - MG Coiled PWR Cord-LH
PCMR	9A669001	Plastic Bin - MG Coiled PWR Cord-RH
RBWS	9A650001	Rear Basket
SBCB	9A652001	Single Basket-Back
SBCF	9A652001	Single Basket-Front
SHEP	3-004-6030	Spiro Holder Equipment Pole Mount
SHLH	9A654001	Scale Holder-LH
SHRH	9A654001	Scale Holder-RH
TBDH	9A661001	Tablet Display Holder
USCB	9A656001	Adjustable Utility Shelf
VESA	9A660001	VESA Mounting Plate
VIPL	9A649001	VOIP Shelf-LH Mount
VIPR	9A649002	VOIP Shelf-RH Mount



## Midmark 6281 Workstation

- 22" height adjustment
- 19.5" W x 14" D worksurface
- Slide-out mouse surface
- Up to 45" reach off wall
- Weight capacity: 0-8 lb
- Folds to 9" flat to wall
- Unit weight: 33 lb
- Shipping weight: 40 lb

### Model

6281-001

### Description

Laptop Wall Mounted Workstation

### Accessories

Plant Installed	Field Installed	
n/a	9A588001	20" Extension Arm
n/a	9A591001	Track Mount Bar Code Holder
n/a	9A590001	Tamper Proof PC Holder
n/a	9A589001	Adjustable PC Holder



## Midmark 6282 Workstation

- 22" worksurface height adjustment
- 25" W x 10" D worksurface
- Max monitor capacity is 20 lb
- 9.5" independent monitor height adjustment
- Independent 180° monitor rotation and 0-45° tilt
- VESA mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight capacity: 5 - 20 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

**Model**  
6282-001

**Description**  
Standard Duty Wall Mounted Workstation

### Accessories

Plant Installed	Field Installed	
n/a	9A593001	Dual Monitor Arm
n/a	9A588001	20" Extension Arm
n/a	9A592001	VESA Mount Bar Code Holder
n/a	9A591001	Track Mount Bar Code Holder
n/a	9A590001	Tamper Proof PC Holder
n/a	9A589001	Adjustable PC Holder



## Midmark 6283 Workstation

- 22" worksurface height adjustment
- 25" W x 10" D worksurface
- Max monitor capacity is 30 lb
- 9.5" independent monitor height adjustment
- Independent 180° monitor rotation and 0-45° tilt
- VESA mount 75 mm and 100 mm
- Folds to 9.47" flat to wall
- Weight capacity: 20-35 lb
- Product weight: 43 lb
- Shipping weight: 50 lb

### Model

6283-001

### Description

Heavy Duty Wall Mounted Workstation

### Accessories

Plant Installed	Field Installed	
n/a	9A593001	Dual Monitor Arm
n/a	9A588001	20" Extension Arm
n/a	9A592001	VESA Mount Bar Code Holder
n/a	9A591001	Track Mount Bar Code Holder
n/a	9A590001	Tamper Proof PC Holder
n/a	9A589001	Adjustable PC Holder

## Synthesis® Cabinetry Collection by Midmark: Limited Warranty

### SCOPE OF WARRANTY

Midmark Corporation warrants to the original purchaser its Midmark Modular cabinetry products (Synthesis® cabinetry Collection) manufactured by Midmark to be free from defects in material and workmanship under normal interior use and service. Midmark's obligation under this warranty is limited to the repair or replacement, at Midmark's option, of the parts or the products the defects of which are reported to Midmark within the applicable warranty period and which upon examination by Midmark prove to be defective.

### APPLICABLE WARRANTY PERIOD

The applicable warranty period, measured from the date of delivery to the original user, shall be as follows:

- A. Five years for all Midmark brand casework products and components (except for the shorter periods in B and C), including door and drawer fronts and casters
- B. Three years for task lights
- C. One year for countertop surfaces and accessories

### EXCLUSIONS

This warranty does not cover sinks, faucets and plumbing accessories. Any warranties on these items are extended directly by the manufacturer of these items to the original purchaser. Information on these manufacturers' warranties will be enclosed with the applicable products. Also, Midmark will furnish copies of any warranties extended by any such manufacturers upon request.

This warranty does not cover and Midmark shall not be liable for: (1) repairs and replacements required because of misuse, abuse, negligence, alteration, accident, freight damage or tampering; (2) matching of color, grain or texture except to commercially acceptable standards; (3) changes in color caused by natural or artificial light; (4) products which are not installed, used or properly cleaned as required in Midmark's written Installation Manual; (5) specially manufactured products; (6) products considered to be of a consumable nature such as light bulbs and surge suppression products; (7) accessories or parts not manufactured by Midmark; (8) charges by anyone (including Midmark's authorized dealers) for adjustments, repairs, replacement parts, installation or other work performed upon or in connection with such products which is not expressly authorized in writing in advance by Midmark.

### EXCLUSIVE REMEDY

Midmark's only obligation under this warranty is the repair or replacement of defective parts. Midmark shall not be liable for any direct, special, indirect, incidental, exemplary or consequential damages or delay.

### NO AUTHORIZATION


No person or firm is authorized to create for Midmark any other obligation or liability in connection with the products.

THIS WARRANTY IS MIDMARK'S ONLY WARRANTY AND IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED. MIDMARK MAKES NO IMPLIED WARRANTIES OF ANY KIND INCLUDING ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. THIS WARRANTY IS LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE PARTS.

\*Warranty information included in product manuals


# Midmark Synthesis® Cabinetry Color-Specific Part Number Explanation and Standard Finish Options

E1RX . 803 . 362



Example: E1 36" H Exam Room Package, Right Hand Sink, No Options, Henna, Concrete Stone

N2A . 803



Example: Dressing Nook, LH Door, Keyed Locks, Henna



Flare Handle  
Serenity Panel

## Frame / Panel Color:

### Solids:

216 Pebble  
490 Fawn  
649 Whisper  
791 Frost

### Metallics:

804 Venus Silver (M)  
805 Pearl Essence (M)  
806 Earthen Bronze (M)  
807 Flat Iron (M)

### Others:

793 Pebble PVC-Free

### Woodgrains:

459 Acorn  
487 Earth  
794 Dune  
795 Flax  
796 Nest  
797 Hewn  
798 Path  
799 Sediment  
800 Radiance  
801 Timber  
802 Storm  
803 Henna

## Laminate Countertop Color:

362 Concrete Stone  
443 Kalahari Topaz  
364 Graphite Nebula  
360 Evening Tigris  
453 Milano Quartz  
290 Green Tigris  
354 Bronze Legacy  
274 Pebble Grey  
711 Colorado Slate



# Exam Room Cabinetry Package E1

Exam Base: 36" H (39" to top of backsplash) x 48" W x 18" D  
(2) Doors w/Waste Drop, (2) 4" Drawers and (2) Doors

Exam Overhead: 48" W x 24" H x 13" D w/3 position Glove,  
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in sink
- Gooseneck facet w/ 4" wrist-blade handles
- Waste can
- Serenity panel style w/flare handles
- Shipping Weight: 327 lb

## Product Information

### Step 1: Configuration

**E1**

### Step 2: Sink Location

**L** Left Hand

**R** Right Hand

### Step 3: Option Package

**X** None

**A** Keyed Lock

### Step 4: Cabinetry Finish\*

**XXX** Solid Color Cabinet Finish  
(Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free  
Cabinet Finish

### Step 5: Countertop Color\*

**XXX** Countertop Color

**Configuration Example:** E1LX.803.362

\*See page 71 for Cabinetry Finish and Countertop Color





# Exam Room Cabinetry Package E2

Exam Base: 36" H (39" to top of backsplash) x 48" W x 18" D  
(2) Doors w/Waste Drop, (3) 4" Drawers and (1) 8" Drawer

Exam Overhead: 48" W x 24" H x 13" D w/3 position Glove,  
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in sink
- Gooseneck facet w/ 4" wrist-blade handles
- Waste can
- Serenity panel style w/flare handles
- Shipping Weight: 331 lb

## Product Information

### Step 1: Configuration

**E2**

### Step 2: Sink Location

**L** Left Hand

**R** Right Hand

### Step 3: Option Package

**X** None

**A** Keyed Lock

### Step 4: Cabinetry Finish\*

**XXX** Solid Color Cabinet Finish  
(Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free  
Cabinet Finish

### Step 5: Countertop Color\*

**XXX** Countertop Color

**Configuration Example:** E2LX.803.362

\*See page 71 for Cabinetry Finish and Countertop Color



# Exam Room Cabinetry Package E3

Overall Exam Base Dimensions: 36" H x 54" W x 18" D w/Support Leg

Exam Base Cabinet: 36" H (39" to top of backsplash) x 24" W x 18" D  
(2) Doors w/Waste Drop

Exam Cart: 33" H x 24" W x 18" D (3) 4" Drawers, (1) 8" Drawer  
w/Kydex Contour top and 60mm casters

Exam Overhead: 24" H x 48" W x 13" D w/(3) position Glove,  
Towel & Cup Dispenser

- 55" W x 19" D Postform laminate top
- Stainless steel drop-in sink
- Gooseneck facet w/ 4" wrist-blade handles
- Waste can
- Serenity panel style w/flare handles
- Shipping Weight: 437 lb

## Product Information

### Step 1: Configuration

**E3**

### Step 2: Sink Location

**L** Left Hand

**R** Right Hand

### Step 3: Option Package

**X** None

**A** Keyed Lock

### Step 4: Cabinetry Finish\*

**XXX** Solid Color Cabinet Finish  
(Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free  
Cabinet Finish

### Step 5: Countertop Color\*

**XXX** Countertop Color

## Configuration Example: E3LX.803.362

\*See page 71 for Cabinetry Finish and Countertop Color



# Exam Room Cabinetry Package E4

Overall Exam Base Dimensions: 36" H (39" to top of backsplash)  
x 48" W x 18" D

Exam Base Cabinet: 36" H x 24" W x 18" D (2) 4" Drawers & (2) Doors  
Wall Mounted Sink: 24" W X 18" D

Exam Overhead: 24" H & 48" W x 13" D w/(3) position Glove,  
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 307 lb

## Product Information

### Step 1: Configuration

**E4**

### Step 2: Sink Location

**L** Left Hand

**R** Right Hand

### Step 3: Option Package

**X** None

**A** Keyed Lock

### Step 4: Cabinetry Finish\*

**XXX** Solid Color Cabinet Finish  
(Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free  
Cabinet Finish

### Step 5: Countertop Color\*

**XXX** Countertop Color

## Configuration Example: E4LX.791.362

\*See page 71 for Cabinetry Finish and Countertop Color



# Exam Room Cabinetry Package E5

Overall Exam Base Dimensions: 36" H (39" to top of backsplash)  
x 54" W x 18" D w/Support Leg

Exam Cart: 33" H x 24" W x 18" D Cart (3) 4" Drawers, (1) 8" Drawer  
w/Kydex Contour Top and 60mm Casters  
Wall Mounted Sink: 24" W X 18" D

Exam Overhead: 24" H & 48" W x 13" D w/(3) position Glove,  
Towel & Cup Dispenser

- 55" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 401 lb

## Product Information

### Step 1: Configuration

**E5**

### Step 2: Sink Location

**L** Left Hand

**R** Right Hand

### Step 3: Option Package

**X** None

**A** Keyed Lock

### Step 4: Cabinetry Finish\*

**XXX** Solid Color Cabinet Finish  
(Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free  
Cabinet Finish

### Step 5: Countertop Color\*

**XXX** Countertop Color

**Configuration Example:** E5LX.803.362

\*See page 71 for Cabinetry Finish and Countertop Color



# Exam Room ADA Cabinetry Package E6

Overall Exam Base Dimensions: 34" H ADA (37" to top of backsplash)  
x 48" W x 18" D

Exam Base Cabinet: 34" H ADA x 18" W x 18" D (2) 4" Drawers  
& (1) 8" Drawer

ADA Wall Mounted Sink: 30" W X 18" D

Exam Overhead: 24" H x 48" W x 13" D w/(3) position Glove,  
Towel & Cup Dispenser

- 49" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 294 lb

\*Customer is responsible for ensuring product meets Federal, state and local code.

## Product Information

### Step 1: Configuration

**E6**

### Step 2: Sink Location

**L** Left Hand

**R** Right Hand

### Step 3: Option Package

**X** None

**A** Keyed Lock

### Step 4: Cabinetry Finish\*

**XXX** Solid Color Cabinet Finish  
(Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free  
Cabinet Finish

### Step 5: Countertop Color\*

**XXX** Countertop Color

**Configuration Example:** E6LX.803.362

\*See page 71 for Cabinetry Finish and Countertop Color



# Exam Room ADA Cabinetry Package E7

Overall Exam Base Dimensions: 34" H ADA (37" to top of backsplash)  
x 60" W x 18" D w/Support Leg

Exam Cart: 31" H x 24" W x 18" D Cart (2) 4" Drawers, (1) 8" Drawer  
w/Kydex® Contour Top and 60 mm Casters  
Wall Mounted Sink: 30" W X 18" D

Exam Overhead: 24" H x 48" W x 13" D w/(3) position Glove,  
Towel & Cup Dispenser

- 61" W x 19" D Postform laminate top
- Stainless steel drop-in wall mounted sink
- Gooseneck facet w/ 4" wrist-blade handles
- Serenity panel style w/flare handles
- Shipping Weight: 294 lb

\*Customer is responsible for ensuring product meets Federal, state and local code.

## Product Information

### Step 1: Configuration

**E7**

### Step 2: Sink Location

**L** Left Hand

**R** Right Hand

### Step 3: Option Package

**X** None

**A** Keyed Lock

### Step 4: Cabinetry Finish\*

**XXX** Solid Color Cabinet Finish  
(Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free  
Cabinet Finish

### Step 5: Countertop Color\*

**XXX** Countertop Color

**Configuration Example:** E7LX.803.362

\*See page 71 for Cabinetry Finish and Countertop Color



\* Shown with Kydex®  
contour top

# Mobile Treatment Cabinet M2

33" H x 27" W x 18" D (2) 4" Drawers & (2) Doors

- Serenity panel style w/flare handles
- Shipping Weight: 155 lb

## Product Information

### Step 1: Configuration

#### M2

### Step 2. Depth

1 18" Deep

### Step 3. Option Package

- A 60mm Casters
- B 4" Quiet Casters
- C 60mm Casters & Keyed Locks
- D 4" Quiet Casters & Keyed Locks

### Step 4. Cabinetry Finish\*

- XXX Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

### Step 5. Countertop Type

- KCN Kydex Contour Top
- SEB Soft Edge Bumper Top
- SEH Soft Edge Handle Top (29"W)

**Configuration Example:** M21A.791.KCN

\*See page 71 for Cabinetry Finish



\* Shown with soft edge bumper top

# Mobile Treatment Cabinet M4

33" H x 27" W x 18" D (3) 4" Drawers & (1) 8" Drawer

- Serenity panel style w/flare handles
- Shipping Weight: 152 lb

## Product Information

### Step 1: Configuration

#### M4

### Step 2: Depth

1 18" Deep

### Step 3: Option Package

- A 60mm Casters
- B 4" Quiet Casters
- C 60mm Casters & Keyed Locks
- D 4" Quiet Casters & Keyed Locks

### Step 4: Cabinetry Finish\*

- XXX Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

### Step 5: Countertop Type

- KCN Kydex® Contour Top
- SEB Soft Edge Bumper Top
- SEH Soft Edge Handle Top (29"W)

**Configuration Example:** M41A.791.SEB

\*See page 71 for Cabinetry Finish





# Mobile Treatment Cabinet M5

33" H x 27" W x 18" D (5) 4" Drawers

- Serenity panel style w/flare handles
- Shipping Weight: 159 lb

\* Shown with soft edge handle top

## Product Information

### Step 1: Configuration

#### M5

### Step 2: Depth

1 18" Deep

### Step 3: Option Package

- A 60mm Casters
- B 4" Quiet Casters
- C 60mm Casters & Keyed Locks
- D 4" Quiet Casters & Keyed Locks

### Step 4: Cabinetry Finish\*

- XXX Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX Woodgrain, Metallic & PVC-Free Cabinet Finish

### Step 5: Countertop Type

- KCN Kydex® Contour Top
- SEB Soft Edge Bumper Top
- SEH Soft Edge Handle Top (29"W)

**Configuration Example:** M51A.791.SEH

\*See page 71 for Cabinetry Finish

# Glove, Towel + Cup Dispenser



\* D2 shown

- 18" W x 8" D x 24" H
- Shipping weight: 37 lb
- Serenity panel style w/flare handle



\* D4 shown

## Product Information

### Step 1. Configuration

- D1** Glove/Towel/Cup Dispenser - (2) Opening RH Door
- D2** Glove/Towel/Cup Dispenser - (2) Opening LH Door
- D3** Glove/Towel/Cup Dispenser - (3) Opening RH Door
- D4** Glove/Towel/Cup Dispenser - (3) Opening LH Door

### Step 2. Option Package

- X** No Lock

### Step 3. Cabinetry Finish\*

- XXX** Solid Color Cabinet Finish (Excludes PVC-Free)
- XXX** Woodgrain, Metallic & PVC-Free Cabinet Finish

### Configuration Example: D2X.803

\*See page 71 for Cabinetry Finish



# Dressing Nook

- Left-hand door pictured
- 24" W x 7" D x 63" H
- Interior mirror, garment hooks and shoe storage included
- Serenity panel style w/flare handle
- Shipping weight: 138 lb



## Product Information

### Step 1. Configuration

**N1** Dressing Nook - RH Door

**N2** Dressing Nook - LH Door

### Step 2. Option Package

**X** No Lock

**A** Keyed Locks

### Step 3. CabinetryFinish\*

**XXX** Solid Color Cabinet Finish (Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free Cabinet Finish

**Configuration Example:** N2X.803

\*See page 71 for Cabinetry Finish



# Waste Receptacle

- 12" W x 12" D x 24" H
- Capacity: 33 liquid qts
- Serenity panel style
- Shipping weight: 35 lb

## Product Information

Step 1. Configuration

**W1** Waste Receptacle

Step 2. Option Package

**X** No Lock

Step 3. Cabinetry Finish\*

**XXX** Solid Color Cabinet Finish (Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free Cabinet Finish

**Configuration Example:** W1X.803

\*See page 71 for Cabinetry Finish



# Instrument Panel

- Instruments not included
- Shipping weight: 17 lb

## Product Information

### Step 1. Configuration

**I1** Instrument Panel - 30" x 20"

**I2** Instrument Panel - 34" x 20"

### Step 2. Option Package

**X** No Lock

### Step 3. Cabinetry Finish\*

**XXX** Solid Color Cabinet Finish (Excludes PVC-Free)

**XXX** Woodgrain, Metallic & PVC-Free Cabinet Finish

**Configuration Example:** I2X.803

\*See page 71 for Cabinetry Finish



# Ritter Waste Receptacle

- With bag retaining ring
- 14" W x 16" D x 25" H
- Shipping weight: 17 lb
- Capacity: 43 liquid qts.
- Color: Pebble

## Model

261-001

## Description

Waste receptacle



## Ritter M11™ Steam Sterilizer

- Automatic Door
- Chamber size: 11" x 18"
- Meets UL 61010-1, 2nd Edition, IEC 61010-2-040, 1st Edition
- Meets CAN/CSA C22.2 No. 61010-1, 2nd Edition
- Meets CSA C22.2 No. 61010-2-040-07 Part 2-040, 1st Edition
- ASME Certified
- Sterilizer weight: 99 lb (44.9 kg)
- Shipping weight: 131 lb (59.4 kg)

### Model

M11-022

### Description

Ritter M11 - large table-top automatic sterilizer with automatic door



## Ritter M9™ Steam Sterilizer

- Automatic Door
- Chamber size: 9" x 15"
- Meets UL 61010-1, 2nd Edition, IEC 61010-2-040, 1st Edition
- Meets CAN/CSA C22.2 No. 61010-1, 2nd Edition
- Meets CSA C22.2 No. 61010-2-040-07 Part 2-040, 1st Edition
- ASME Certified
- Sterilizer weight: 73 lb (33.1 kg)
- Shipping weight: 81 lb (36.7 kg)

### Model

M9-022

### Description

Ritter M9 - fully automatic sterilizer with automatic door



## Ritter M9D™

### Manual Door

- Chamber size: 9" x 15"
- Meets UL 61010-1, 2nd Edition, IEC 61010-2-040, 1st Edition
- Meets CAN/CSA C22.2 No. 61010-1, 2nd Edition
- Meets CSA C22.2 No. 61010-2-040-07 Part 2-040 1st Edition
- ASME Certified
- Sterilizer weight: 73 lb (33.1 kg)
- Shipping weight: 81 lb (36.7 kg)

#### Model

M9D-022

#### Description

Ritter M9D - automatic sterilizer with manual door



## Ritter M3™

### Steam Sterilizer

- Automatic Door
- Chamber size: 12.1" L x 7.6" W x 1.6" D (0.49 gal.)
- Meets UL 61010-1 2nd edition, IEC 61010-2-040
- Meets CAN/CSA C22.2 No. 60110-1 2nd Edition
- Meets CSA C22.2 No. 61010-2-040
- FCC Part 15, Sub-part B
- ASME Certified
- Sterilizer weight: 71 lb (32.2 kg)
- Shipping weight: 80 lb (36.3 kg)

#### Model

M3-003

#### Description

Ritter M3 - compact, fast, automatic sterilizer



## Sterilizer Accessories

Model	Description
002-0396-00	Speed-Clean (16 oz. single bottle) (C)
002-0396-05	Speed-Clean (12 - 16 oz. bottles)
9A226001	M9/M11 Pouch Rack (D)
9A259001	M9/M11 Impact Dot Matrix Printer (B)
9A307001	M9/M11 Cool Hand Tool (Tray and Cassette Removal Tool) (A)
9A401001	M3 Printer (E)
9A402001	M3 Door Tray (F)
9A404001	M3 Top Cover Protector (G)

(A) Cool Handtool

8.5" L



E) M3 Printer



(B) Impact Dot Matrix Printer

7.5" L x 4.4" W x 2.25" H



(F) M3 Door Tray



(C) SpeedClean


8" H



(G) M3 Top Cover Protector



(D) Pouch Rack





## Single VistaCool™

### Direct-to-Drain System

- Designed to be compatible with the Ritter M3 sterilizer
- Tank dimensions w/fittings: 7.6" W x 7.6" L x 13.75" H  
(19.3 cm x 19.3 cm x 34.9 cm)
- Tank volume: 1.50 gal (5.70 liters)
- Tank weight w/water and fittings: 15.0 lb (6.8 kg)
- Approximate shipping weight: 7.0 lb (3.2 kg)
- Max coolant water pressure: 100 psi (6.9 bar)

#### Model

9A586001

#### Description

VistaCool single unit for connecting to one M3 sterilizer



## Double VistaCool™

### Direct-to-Drain System

- Designed to be compatible with Ritter M11 or M9 sterilizer
- Tank dimensions w/fittings: 4.0" W x 13.0" L x 17.0" H  
(10.2 cm x 33.1 cm x 43.2 cm)
- Tank volume: 3.00 gal (11.40 liters)
- Tank weight w/water and fittings: 28.0 lb (12.7 kg)
- Approximate shipping weight: 10.0 lb (4.5 kg)
- Max coolant water pressure: 100 psi (6.9 bar)

#### Model

9A586002\*

#### Description

VistaCool double unit for connecting to any two current Ritter sterilizers

\*Required adapter kits included. Ritter M9 and M11 only compatible with the double unit.

## Accessories

#### Model

#### Description

9A586004

Air Gap Kit (needed only if required by local code)

9A586005

Single VistaCool adjustable tank stand

9A586006

Double VistaCool adjustable tank stand

VistaCool™ is a trademark of Vista Research Group, Ashland, OH



## QuickClean® Ultrasonic Cleaner QC6/QC6R

- Tank Capacity: 6.6 gal/25L
- Tank Size: 19.9" L x 11.8" W x 7.9" H (50.5 x 30 x 20 cm)
- Overall Size: 21.6" L x 16.5" W x 14.3" H (55 x 42 x 36.5 cm)
- Table Top Unit Weight: 29.6 lb (13.44 kg)
- Table Top Unit Shipping Weight: 32.6 lb (14.78 kg)
- Recessed Unit Weight: 32.9 lb (14.92 kg)
- Recessed Unit Shipping Weight: 35.8 lb (16.26 kg)
- Safety basket INCLUDED
- UL 61010-1:2012
- CAN/CSA C22.2 No. 61010-1:2012

Model	Description
QC6-01	QuickClean Ultrasonic Cleaner – 6.6 gal Tabletop Unit
QC6R-01	QuickClean Ultrasonic Cleaner – 6.6 gal Recessed Unit



## QuickClean® Ultrasonic Cleaner QC3/QC3R

- Tank Capacity: 3.3 gal/12.5L
- Tank Size: 13" L x 11.8" W x 5.9" H (33 x 30 x 15 cm)
- Overall Size: 14.7" L x 16.5" W x 10.4" H (37.5 x 42 x 26.5 cm)
- Table Top Unit Weight: 18.6 lb (8.42 kg)
- Table Top Unit Shipping Weight: 20.9 lb (9.50 kg)
- Recessed Unit Weight: 22.3 lb (10.10 kg)
- Recessed Unit Shipping Weight: 24.6 lb (11.18 kg)
- Safety basket INCLUDED
- UL 61010-1:2012
- CAN/CSA C22.2 No. 61010-1:2012

Model	Description
QC3-01	QuickClean Ultrasonic Cleaner – 3.3 gal Tabletop Unit
QC3R-01	QuickClean Ultrasonic Cleaner – 3.3 gal Recessed Unit



## QuickClean® Ultrasonic Cleaner QC1

- Tank Capacity: 1.2 gal/4.5L
- Tank Size: 11.8" L x 5.9" W x 5.9" H (30 x 15 x 15 cm)
- Overall Size: 13.5" L x 10.2" W x 10.4" H (30 x 26 x 26.5 cm)
- Table Top Unit Weight: 12.7 lb (5.76 kg)
- Table Top Unit Shipping Weight: 14.6 lb (6.64 kg)
- Safety basket INCLUDED
- UL 61010-1:2012
- CAN/CSA C22.2 No. 61010-1:2012

Model	Description
QC1-01	QuickClean Ultrasonic Cleaner – 1.2 gal Tabletop Unit

---

## Ultrasonic Parts and Accessories

9A612001	QC1 QuickClean Beaker Accessory Kit
9A613001	QC3/QC3R QuickClean Beaker Accessory Kit
9A614001	QC6/QC6R QuickClean Beaker Accessory Kit
9A298001	Midmark Enzymatic Cleaner (32 oz)
9A296001	Midmark General Purpose Cleaner (32 oz)
9A297001	Midmark Tartar and Stain Remover (32 oz)
002-10007-00	QC1 Safety Basket (1 basket comes standard with a QC1 unit)
002-10008-00	QC3/QC3R Safety Basket (1 basket comes standard with QC3/QC3R unit)
002-10009-00	QC6/QC6R Safety Basket (1 basket comes standard with QC6/QC6R unit)



Designing better care.™

**CARB 93120.2 Phase 2 Compliant**

**Midmark is an ISO 13485 and ISO 9001 Certified Company.**  
For more information or a demonstration, contact your Midmark  
dealer or call: 1-800-MIDMARK Fax: 1-800-365-8631  
or visit our website at [midmark.com](http://midmark.com)

© 2019 Midmark Corporation  
Midmark Corporation, Miamisburg, Ohio USA  
Products subject to improvement changes without notice  
Litho in U.S.A. MKT00176 (2/19)

  
**midmark**™